

# C2 Camera System



The camera and light source combination

# C2 CAMERA & LIGHT SOURCE COMBINATION

This integrated FullHD-Camera and LED Light Source is a combined 2-in-1 system. The concept is specially designed for the use in the outpatient market.

The C2 support rigid scopes and fiber-endoscopes and with the small footprint it provides an economical and flexible solution for endoscopic diagnostic and surgical interventions in Gynecology, Urology or ENT.

## Features:

- Compact and integrated 2-in-1 design
- Versatile use for rigid scopes and fiberoptics
- LED lamp – exceptionally long lifespan
- Low noise ventilation system
- Easy storage and documentation of videos and pictures
- Optional TV couplers with fixed focal length or with zoom available

**FULL HD**



Integrated recording functions



Rigid and fiber endoscopes supported



Easy to use



250-015-030  
Zoom TV Adapter

## DELIVERY CONTAINS:

- Main Unit
- Camera head
- HDMI cable
- Mains cable
- Operating manual

C2 Camera and Light Source Combination

**300-002-200**

## Technical data

Line voltage	100-240 VAC
Power consumption	60VA
Image sensor	1/2" MOS
Resolution	1920 x 1080 (FullHD)
Index ratio	16:9
Lamp	LED (50W comparable)
Video outputs	1x HDMI 1x DVI 1x SDI 1x VIDEO 1x USB (A) 1x Remote IN (Footswitch) 1x RJ45 (Service)
Dimensions	Control unit: 295 x 128 x 279mm (WxHxD) Weight control unit: approx. 3kg Weight camera head: 180g

TV Coupler, Endoscope and Light Guide cable not included and can be ordered separately.  
NOTE: Endoscopes smaller than Ø 2.9mm require 250-015-030 Zoom TV Adapter.



RZ Medizintechnik GmbH

Unter Hasslen 20/22  
78532 Tuttlingen, Germany

Phone: +49 7462 9470-0

[sales@rz-medizintechnik.com](mailto:sales@rz-medizintechnik.com)  
[www.rz-medizintechnik.com](http://www.rz-medizintechnik.com)

