

360° Rotatable,
color-coded
Instruments




360° ROTATABLE INSTRUMENTS

Dedicated Instruments for Urology and Hysteroscopy

The unique design and features of the rotating color coded instruments have been designed in cooperation with Dr Accardi. Multiple years of experience in using dedicated instruments have led to this design and functionality.

The automatic holding mechanism combined with the 360° rotation of the jaw and the light design will improve your comfort.



Color-coded wheel to indicate the different jaws

Overload protection helps to prevent damage

Surface lotus effect increase corrosion resistance

360° rotation just with your finger tip



Ergonomic, lightweight handle







Available in 4 and 5 Charr. diameter - rigid, semi rigid and flexible

Automatic closing mechanism - active holding of tissue due to internal spring system

360° ROTATABLE INSTRUMENTS

Instruments for Urology

Flexible and rigid instruments of 600mm length and different designs with various jaws for many indications in Urology.




Knob Color		5 Charr.	4 Charr.	
■	rigid	253-950-005	253-950-004	Stone Grasping Forceps, double action jaws
	flexible	253-951-005	253-951-004	
■	rigid	253-952-005	253-952-004	Stone Grasping Forceps, large stones, double action jaws
	flexible	253-953-005	253-953-004	
■	rigid	253-954-005	253-954-004	HAMADA Long Stone Grasping Forceps, double action jaws
	flexible	253-955-005	253-955-004	
■	rigid	253-956-005	253-956-004	Biopsy Forceps, double action jaws
	flexible	253-957-005	253-957-004	
■	rigid	253-958-005	253-958-004	Scissors, single action jaws
	flexible	253-959-005	253-959-004	



360° ROTATABLE INSTRUMENTS

Instruments for Hysteroscopy

Semi-rigid instruments with dedicated jaws for your surgical hysteroscopy. The color codes will help always having the perfect instruments for your clinical needs.



	Knob Color	400mm	340mm	
		351-116-305	351-116-335	Grasping Forceps with cup, serrated, double action
		351-118-305	351-118-335	Scissors, sharp tips, single action
		351-100-305	351-100-335	Scissors, blunt tips, single action
		351-126-305	351-126-335	Biopsy Punch Forceps, single action
		351-120-305	351-120-335	Biopsy Forceps, double action
		351-117-305	351-117-335	Tenaculum Forceps, 1x1 teeth, double action
		351-115-305	351-115-335	SHARK Tenaculum Grasping Forceps
		351-119-305	351-119-335	Accardi HARPOON Grasping Forceps, with barbs





RZ Medizintechnik GmbH

Unter Hasslen 20
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com



Medica Code Number 021542