




ArthVANCE<sup>®</sup>  
SHOULDER  
INSTRUMENTS

The 'ArthVANCE' logo features a stylized 'A' with a teal brushstroke effect. Below it, the text 'ArthVANCE' is in a grey, sans-serif font with a registered trademark symbol. Underneath, the words 'SHOULDER' and 'INSTRUMENTS' are stacked in a bold, black, all-caps, sans-serif font.

# ArthVANCE® SHOULDER INSTRUMENTS


## SUTURE PASSER

Simultaneously one-handed passing and retrieving of suture


- Efficient tissue bite - wide jaw opening and deep bite
- Minimal trauma at tissue penetration
- Reusable instrument and reliable construction




Automatic trap loop hold




Strong grip with 4-finger trigger



Slim design



Suture Passer with trap door, frontal opening, without ratchet, 170mm working length  
225-820-217



Suture Passer with trap door, frontal opening, with ratchet, 170mm working length  
225-820-213




Needle with lateral opening, 0.3 x 240mm  
225-820-219


## SUTURE CUTTER / KNOT PUSHER



Knot Pusher, single-hole, for sliding knots and half-stitches, modified distal tip for easy slipknots, 175mm  
225-259-000




Suture Cutter, axial handle, Ø 4.2 mm, v-notch, left opened, 130 mm working length  
225-820-100



Suture Cutter, ring handle, Ø 4.2 mm, v-notch, left opened, 130 mm working length  
225-820-101

# ArthVANCE® SHOULDER INSTRUMENTS



Fiber Retriever, straight,  $\varnothing$  1.5mm pin  
155mm working length  
225-820-259



Debridement Rasp, razor sharp,  
ergonomic PEEK handle  
225-600-630



Tissue Elevator, precise, 15° angled,  
ergonomic PEEK handle  
225-600-650

# SUTURE INSTRUMENTS

For rotator and labral repair



## Smart Hook



Left curved jaw

**A** 225-820-307

**B** 225-840-307



Right curved jaw

**A** 225-820-309

**C** 225-840-309

## Sixer



Twisted left

**A** 225-820-306

**B** 225-840-306




Twisted right

**A** 225-820-308

**C** 225-840-308

# ArthVANCE® SHOULDER INSTRUMENTS

## Birdy



Straight jaw, upcurved shaft

**A** 225-820-400

**B** 225-840-400



20° upwards angled

**A** 225-820-402

**B** 225-840-402



40° upwards angled

**A** 225-820-404

**B** 225-840-404



45° right angled

**A** 225-820-406

**C** 225-840-406



45° left angled



**A** 225-820-408

**B** 225-840-408



Rotator Cuff Simulator (available on request)

225-225-225



Please do not test the Suture Passer without adequate material in between the jaw.  
The needle may break due to the friction.



## RZ Medizintechnik GmbH

Unter Hasslen 20  
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50

[sales@rz-medizintechnik.com](mailto:sales@rz-medizintechnik.com)  
[www.rz-medizintechnik.com](http://www.rz-medizintechnik.com)

