

INDEX

Sets	2 - 20
Cystoscopy	21 - 38
Resectoscopy	39 - 64
Laser Enucleation	65 - 77
Urethrotome	78 - 82
Lithotripsy	83 - 87
Nephroscopy	88 - 96
Uretero - Renoscopy	97 - 102
Pediatric	103 - 112
Devices and Flexible Instruments	113 - 115
Accessories	116 - 121
Cleaning	122 - 130
Equipment	131 - 146




Products marked with this symbol are not available with RZ-marking, and not the CE sign/CE number of RZ.



RZ-Topseller products

All products marked with this green star are available on stock.

CONFIGURATIONS by RZ Cystoscope Set for adults




253-900-900 Cystoscope Set for adults

1	RZ Cystoscope, 30°, Ø 4.0mm	253-804-030	1 pc.
2	Cysto-Urethroscope Sheath 17 Charr.	253-300-017	1 pc.
3	Cysto-Urethroscope Sheath 22 Charr.	253-300-022	1 pc.
4	Bridge with 2 channels for 1 instrument with max. 12 Charr. or 2 instruments with max. 8 Charr.	253-000-102	1 pc.
5	Albarran deflecting mechanism with 2 working channels and ratchet	253-000-105	1 pc.
6	Flexible Grasping Forceps, serrated jaws, double action	351-114-107	1 pc.
7	Biopsy Forceps, oval cup, serrated jaws, double action	351-120-107	1 pc.
8	Scissors, single action	351-100-105	1 pc.
9	Spoon Biopsy Forceps with round cup, double action	253-420-100	1 pc.
10	Dormia Metal Stone Basket, flexible, 4 wires, 600mm working length, 5 Charr.	253-585-005	1 pc.
11	Ball Electrode 530mm working length, 5 Charr.	253-370-005	1 pc.
12	Ball Electrode 530mm working length, 7 Charr.	253-370-007	1 pc.
13	Monopolar Cable with 8 mm plug for Valleylab HF-unit	253-610-003	1 pc.
14	Elik Evacuator with turn lock connection	253-500-000	1 pc.
15	Bladder Syringe, 150cc	253-000-950	1 pc.




Cystoscope Set for children



253-900-918 Cystoscope Set for children

1	RZ Pediatric Fiber Cystoscope, 30°, Ø 2.0mm	253-219-130	1 pc.
2	Pediatric Cystoscope Sheath 7.5 Charr.	253-200-075	1 pc.
3	Pediatric Cystoscope Sheath 9.5 Charr.	253-200-095	1 pc.
4	Alligator grasping forceps, 28cm length, 3 Charr.	253-202-003	1 pc.
5	Biopsy forceps, 28cm length, 3 Charr.	253-204-003	1 pc.
6	Ball electrode, 25cm length, 3 Charr.	253-370-253	1 pc.
7	Needle electrode, 25cm length, 3 Charr.	253-371-253	1 pc.
8	Monopolar Cable with 8 mm plug for Valleylab HF-unit	253-610-003	1 pc.

CONFIGURATIONS by RZ Cystoscope Set for (newborn) babies




253-900-919 Cystoscope Set for (newborn) Babies

1	Operating Cysto-Urethrotome 4.5 Charr., 5° direction of view	253-200-600	1 pc.
2	Stone Grasping Forceps, double action, 3 Charr.	253-203-128	1 pc.
3	Biopsy Forceps, double action, 3 Charr.	253-203-328	1 pc.
4	Scissors, single action, 3 Charr.	253-203-428	1 pc.
5	Ball Electrode, 3 Charr.	253-370-253	1 pc.
6	Needle electrode, 3 Charr.	253-371-253	1 pc.
7	Hook Knife, 3 Charr.	253-200-200	1 pc.
8	Injection Cannula, 3 Charr.	253-200-210	1 pc.



26 Charr. Monopolar Resectoscope Set for adults




253-900-901 26 Charr. Monopolar Resectoscope Set for adults

1	RZ Cystoscope, 30°, Ø 4.0mm	253-804-030	1 pc.
2	Passive Titanium working element	253-000-200	1 pc.
3	Continuous Flow Resectoscope Sheath, 26 Charr., with obturator	253-000-256	1 pc.
4	Visual Obturator for 4mm RZ Cystoscopes	253-000-084	1 pc.
5	Cutting Loop, angled, for 12°/30° telescope	253-300-201	1 pc.
6	Coagulation Electrode, pointed	253-300-203	1 pc.
7	Coagulation Electrode, ball shape, Ø 3mm	253-300-205	1 pc.
8	Coagulation Electrode, ball shape, Ø 5mm	253-300-206	1 pc.
9	Sterilization tubes for electrodes	253-300-999	1 pc.
10	Monopolar Cable with 8 mm plug for Valleylab HF-unit	253-610-003	1 pc.
11	Ellik Evacuator with turn lock connection	253-500-000	1 pc.
12	Bladder Syringe, 150cc	253-000-950	1 pc.



Monopolar Resectoscope Set for children




253-900-905 Monopolar Resectoscope Set for children

1	RZ Pediatric Fiber Cystoscope, 0°, Ø 2.0mm	253-219-100	1 pc.
2	Passive Titanium working element	253-250-200	1 pc.
3	Resectoscope Sheath 11 Charr. with obturator	253-250-211	1 pc.
4	Cutting Loop Electrode, 11 Charr.	253-251-111	1 pc.
5	Coagulation Electrode, 11 Charr.	253-251-112	1 pc.
6	Hook Electrode with ball, 11 Charr.	253-251-113	1 pc.
7	Needle Electrode, pointed, 11 Charr.	253-251-115	1 pc.
8	Monopolar Cable with 8 mm plug for Valleylab HF-unit	253-610-003	1 pc.



26 Charr. Bipolar (HYBRID) Resectoscope Set for adults




253-900-911 26 Charr. Bipolar (HYBRID) Resectoscope Set for adults

1	RZ Cystoscope, 30°, Ø 4mm	253-804-030	1 pc.
2	Passive HYBRID working element	253-000-400	1 pc.
3	Continuous Flow Resectoscope Sheath, 26 Charr.	253-000-256	1 pc.
4	Cutting Electrode	253-400-201	1 pc.
5	Collins Coagulation Electrode	253-400-203	1 pc.
6	Ball Electrode	253-400-205	1 pc.
7	Visual Obturator, with 6 Charr. working channel	253-000-094	1 pc.
8	Sterilization tube for electrodes	253-300-999	1 pc.
9	SDS Cable, 4.5m, bipolar	700-500-453	1 pc.



22 Charr. Bipolar (HYBRID) Resectoscope Set




22 Charr. Hybrid Resectoscope Set

1	RZ Rod Lens Cystoscope, 30°, Ø 2.9mm	351-829-030	1 pc.
2	Passive Titanium Working Element	253-000-410	1 pc.
3	Resectoscope Sheath, 22 Charr., rotating inner tube, for continuous suction / irrigation	253-000-352	1 pc.
4	Loop Electrode, 30° angled	253-000-541	1 pc.
5	Ball Electrode	253-000-545	1 pc.
6	Knife Electrode	253-000-543	1 pc.
7	Storage container for knives and electrodes, metal, perforated	253-300-999	1 pc.
8	SDS Cable, 4.5m, bipolar	700-500-453	1 pc.




18.5 Charr. Bipolar (HYBRID) Resectoscope Set



253-900-921 18.5 Charr. Bipolar (HYBRID) Resectoscope Set for adults and adolescents

1	RZ Cystoscope, 30°, Ø 2.9mm	253-829-030	1 pc.
2	Passive HYBRID working element	351-400-200	1 pc.
3	Continuous Flow Resectoscope Sheath, 18.5 Charr.	351-000-185	1 pc.
4	Loop Electrode	351-000-201	1 pc.
5	Knife Electrode	351-000-203	1 pc.
6	Ball Electrode	351-000-205	1 pc.
7	Visual Obturator, with 5 Charr. working channel	351-000-180	1 pc.
8	Sterilization tube for electrodes	253-300-999	1 pc.
9	SDS Cable, 4.5m, bipolar	700-500-453	1 pc.


CONFIGURATIONS by RZ *slimHoLEP_k* Laser Resectoscope Set according to Dr. Felipe C. A. de Figueiredo



Laser Resectoscope Set according to Dr. Felipe C. A. de Figueiredo

1	RZ Rod Lens Cystoscope, 30°, Ø 2.9mm	351-829-030	1 pc.
2	RZ Laser Working Element, passive, for fibres up to 0.8mm	253-000-319	1 pc.
3	Resectoscope Sheath, 22 Charr., rotating inner tube, for continuous suction / irrigation	253-000-352	1 pc.
4	RZ LL-Connection, rotatable, for use with Laser Probe, short	253-000-301	1 pc.
5	Luer-Lock Tuohy Borst Adapter (Ø 0.6 - 1.4mm)	300-011-184	1 pc.
6	Morcroscope, 22 Charr., 220mm working length, 5mm working channel	253-905-220	1 pc.
7	Sealing Cap for Morcroscope, package of 10 pieces	253-904-220	1 pc.

CONFIGURATIONS by RZ *ultra slim HoLEP_K* Laser Resectoscope Set according to Dr. Felipe C. A. de Figueiredo




Laser Resectoscope Set according to Dr. Felipe C. A. de Figueiredo

1	RZ Rod Lens Cystoscope, 30°, Ø 2.9mm	351-829-030	1 pc.
2	RZ Laser Working Element, passive, for fibres up to 0.8mm, with push button	253-000-318	1 pc.
3	Resectoscope Continuous Flow Sheath, 18.5 Charr.	351-000-185	1 pc.
4	RZ LL-Connection, rotatable, for use with Laser Probe, short	253-000-301	1 pc.
5	Luer-Lock Tuohy Borst Adapter (Ø 0.6 - 1.4mm)	300-011-184	1 pc.
6	Morcscope, 22 Charr., 220mm working length, 5mm working channel	253-905-220	1 pc.
7	Continuous Flow Sheath, 22 Charr.	253-901-220	1 pc.
8	Sealing Cap for Morcscope, package of 10 pieces	253-904-220	1 pc.



HoLEP Set




253-900-920 HoLEP Set

1	RZ Cystoscope, 30°, Ø 4.0mm	253-804-030	1 pc.
2	RZ Laser Working Element, passive, for laser-probes	253-000-312	1 pc.
3	Continuous Flow Resectoscope Sheath, 26 Charr.	253-000-256	1 pc.
4	RZ LL-Connection	253-000-302	1 pc.
5	Luer-Lock Tuohy Borst Adapter	300-011-184	1 pc.
6	Nephroscope 6°	253-900-121	1 pc.
7	Adapter for Resectoscope outer sheath 253-000-256	253-428-100	1 pc.



Urethrotome Set for adults




253-900-902 Urethrotome Set for adults

1	RZ Cystoscope, 0°, Ø 4.0mm	253-804-000	1 pc.
2	Passive Titanium Working Element	253-000-200	1 pc.
3	Urethrotome Sheath oval shape, 21 Charr.	253-500-221	1 pc.
4	Hollow Obturator with channel for guide wires	253-500-321	1 pc.
5	Bridge with channel for 5 Charr. accessories	253-500-325	1 pc.
6	Supplementary Sheath, half round	253-500-230	1 pc.
7	Sachse knife, straight	253-500-561	3 pcs.
8	Half round knife, 180° cutting	253-500-563	2 pcs.
9	WaveCut knife, straight	253-500-564	1 pc.
10	Sterilization tube for electrodes	253-300-999	1 pc.
11	Bladder Syringe, 150cc	253-000-950	1 pc.
12	OTIS Dilation Urethrotome, complete	253-500-100	1 pc.




Urethrotome Set for children



253-900-922 Urethrotome Set for children

1	RZ Pediatric Fiber Cystoscope, 0°, Ø 2.0mm	253-219-100	1 pc.
2	Passive Working Element	253-250-200	1 pc.
3	Pediatric Urethrotome sheath, 10 Charr.	253-265-010	1 pc.
4	Bridge with 1 channel	253-255-001	1 pc.
5	Knife with straight cutting edge	253-265-101	1 pc.
6	Round Knife	253-265-102	1 pc.
7	Hook Knife	253-265-104	1 pc.

CONFIGURATIONS by RZ Uretero-Renoscope Set




253-900-903 Uretero-Renoscope Set

1	Advanced Semi-Rigid Uretero-Renoscope	253-905-540	1 pc.
2	Stone Grasping Forceps	253-950-005	1 pc.
3	Stone Grasping Forceps for large stones	253-952-005	1 pc.
4	HAMADA Long Stone Grasping Forceps	253-954-005	1 pc.
5	Biopsy Forceps	253-956-005	1 pc.
6	Dormia Metal Stone Basket, flexible	253-585-005	1 pc.



Nephroscope Set, Mini PCNL




253-900-909 Nephroscope Set, Mini PCNL

1	Mini Nephroscope, 12°	253-900-140	1 pc.
2	15 Charr. Nephroscope Sheath	253-900-141	1 pc.
3	One-Step Dilator	253-900-146	1 pc.
4	2-piece Sealing Adapter	253-900-145	1 pc.
5	16.5 Charr. Nephroscope Sheath	253-900-142	1 pc.
6	One-Step Dilator	253-900-147	1 pc.
7	21 Charr. Nephroscope Sheath	253-900-143	3 pcs.
8	One-step Dilator	253-900-148	3 pcs.
9	Stone Grasping forceps	253-950-455	1 pc.
10	Stone Grasping Forceps for large stones	253-952-455	1 pc.
11	Biopsy Forceps	253-956-455	1 pc.



Nephroscope Set, Compact PCNL




253-900-923 Nephroscope Set, Compact PCNL

1	Compact Nephroscope	253-900-830	1 pc.
2	Shaft with Obturator	253-900-888	1 pc.
3	Stone Grasping forceps	253-950-455	1 pc.
4	Stone Grasping Forceps for large stones	253-952-455	1 pc.
5	Biopsy Forceps	253-956-455	1 pc.




Nephroscope Set, Standard PCNL



253-900-924 Nephroscope Set, Standard PCNL

1	Nephroscope, 6°	253-900-121	1 pc.
2	Operating Sheath 26 Charr.	253-900-130	1 pc.
3	Hollow Obturator 26 Charr.	253-900-135	1 pc.
4	Amplatz metal Dilator Set	253-900-155	1 pc.
5	Amplatz metal Dilator 27 Charr.	253-900-156	1 pc.
6	Amplatz metal Dilator 30 Charr.	253-900-157	1 pc.
7	Puncture Cannula with inner and outer cannula	253-900-201	3 pcs.
8	Dilation Cannula, Ø 3mm	253-900-202	3 pcs.
9	Stone Grasping Forceps, fenestrated	253-900-250	1 pc.
10	Stone Grasping Forceps, coarse serrated	253-900-252	1 pc.

CONFIGURATIONS by RZ Lithotripsy Set, mechanical



253-900-925 Lithotripsy Set, mechanical

1	RZ Cystoscope, 0°, Ø 4.0mm	253-804-000	1 pc.
2	RZ Cystoscope, 30°, Ø 4.0mm	253-804-030	1 pc.
3	Mauermayer Stone Punch	253-755-010	1 pc.
4	Punch Sheath with Central Valve, 24 Charr.	253-755-020	1 pc.
5	Visual Obturator with 7 Charr. working channel	253-755-030	1 pc.
6	Cysto-Urethroscope Sheath 25 Charr.	253-300-025	1 pc.
7	Stone Crushing Forceps for small stones, single action	253-760-000	1 pc.
8	Ellik Evacuator with turn lock connection	253-500-000	1 pc.



Lithotripsy Set, electro-pneumatic



1



2



3




4 – 9

253-900-926 Lithotripsy Set, electro-pneumatic

1	L 100 Lithotripter	253-750-100	1 pc.
2	Pneumatic Footswitch	253-750-150	1 pc.
3	Hand piece, anodised	253-750-200	1 pc.
4	Litho Probe, 0.6mm	253-850-006	2 pcs.
5	Litho Probe, 0.8mm	253-850-008	2 pcs.
6	Litho Probe, 1.0mm	253-850-010	2 pcs.
7	Litho Probe, 1.2mm	253-850-012	2 pcs.
8	Litho Probe, 1.6mm	253-850-016	2 pcs.
9	Litho Probe, 2.0mm	253-850-020	2 pcs.

CYSTOSCOPY




4K UHD CYSTOSCOPES

- 4K Ultra High Definition resolution
- Ø 4mm, 302mm working length
- Available in 0°, 12°, 30° and 70° direction of view
- Autoclavable
- Made in Germany



RZ 4K Cystoscope, Ø 4mm, 302mm working length, 362mm overall length




FULL HD CYSTOSCOPES


- Colour coded rings identify the direction of view
- Ø 4mm, 302mm working length
- 3-step light post adapters for the most common light cables
- Sapphire glasses to protect lens system
- Optimized glass fibre arrangements for an uniform image brightness
- Fully autoclavable at 134°C and 2.3bar



RZ Cystoscope, Ø 4mm, 302mm working length, 356mm overall length




Sheaths for use with Ø4mm Cystoscopes



Standard Cysto-Urethroscope sheaths with 2 lateral stopcocks and standard obturator with atraumatic tip

- 253-300-017★ Cysto-Urethroscope Sheath 17 Charr. ●
- 253-300-019★ Cysto-Urethroscope Sheath 19 Charr. ●
- 253-300-020★ Cysto-Urethroscope Sheath 20 Charr. ●
- 253-300-022★ Cysto-Urethroscope Sheath 22 Charr. ●
- 253-300-025★ Cysto-Urethroscope Sheath 25 Charr. ○



Obturator for use with Cystoscope sheaths



Standard Obturators with atraumatic tip and locking mechanism for all RZ Cysto-Urethroscope sheaths

253-000-117	Obturator 17 Charr. ●
253-000-119	Obturator 19 Charr. ●
253-000-120	Obturator 21 Charr. ●
253-000-122	Obturator 23 Charr. ●
253-000-125	Obturator 25 Charr. ○




Visual Obturators with locking mechanism for all 4mm RZ Cystoscopes

253-000-217	Obturator 17 Charr. ●
253-000-219	Obturator 19 Charr. ●
253-000-220	Obturator 20 Charr. ●
253-000-222	Obturator 22 Charr. ●
253-000-225	Obturator 25 Charr. ○




253-315-000
Hamada Cystoscope Adapter for use with RZ Cystoscope sheath

Cystoscope Bridge for use with Cysto-Urethroscope sheaths




- 1** 253-000-100* without channel
- 2** 253-000-101* with 1 channel for 1 instrument with max. 12 Charr.
- 3** 253-000-102* with 2 channels for 1 instrument with max. 12 Charr. or 2 instruments with max. 8 Charr.

Albarran deflecting mechanism for use with 19-25 Charr. RZ Cysto-Urethroscope sheaths



- 253-000-105* Albarran deflecting mechanism with 2 working channels **and ratchet**
- 253-000-107 Albarran deflecting mechanism with 1 working channel **and ratchet** (without illustration)




- 253-000-106 Albarran deflecting mechanism with 2 working channels **without ratchet** (without illustration)
- 253-000-108 Albarran deflecting mechanism with 1 working channel **without ratchet**

Sheath - Instrument combinations

SHEATHS	COLOR	ALBARRAN		BRIDGE	
		1 working channel	2 working channels	1 working channel	2 working channels
25 Charr.	○	1 x 10 Charr.	2 x 8 Charr.	1 x 12 Charr.	2 x 8 Charr.
22 Charr.	●	1 x 9 Charr.	2 x 6 Charr.	1 x 10 Charr.	2 x 7 Charr.
20 Charr.	●	1 x 6 Charr.	2 x 5 Charr.	1 x 7 Charr.	2 x 6 Charr.
19 Charr.	●	1 x 5 Charr.	2 x 4 Charr.	1 x 6 Charr.	2 x 5 Charr.
17 Charr.	●	-	-	1 x 5 Charr.	-

Flexible Instruments for use with cystoscopes



Working length: 400mm
Rotatable and detachable



Grasping Forceps, serrated jaws,
double action

351-114-105* 5 Charr.
351-114-107* 7 Charr.
351-114-109 9 Charr.



Biopsy Forceps, oval cup,
double action

351-120-105* 5 Charr.
351-120-107* 7 Charr.
351-120-109 9 Charr.



Serrated Biopsy Forceps, oval cup
with aggressive serration, double action

351-124-107 7 Charr.



Scissors, single action

351-100-105* 5 Charr.
351-100-107 7 Charr.
351-100-109 9 Charr.



DJ Forceps, double action

351-107-107 7 Charr.
351-107-109 9 Charr.



Grasping Forceps, serrated jaw,
double action, 730mm working length

253-958-105 5 Charr.



Biopsy Forceps, double action,
730mm working length

253-957-105 5 Charr.




The next generation instruments in Endo-Urology

- ▼ Detachable – one handle for all inserts
- ▼ Color code – the easy identification of the tip
- ▼ Two handles – self holding and standard
- ▼ Easy connection – the reliable ball tip technique
- ▼ Light weight handle – the perfect fit for your fingers
- ▼ Make cleaning easy – flushing port for the complete instrument
- ▼ 360° rotation - just with your finger tip
- ▼ All parts can be ordered separately



PATENTED
10 2022 107 371.8






Multiple tips design




Color coded inserts



Unique, detachable system



Flushing port for easy cleaning



Two handles - multiple options


FORCEPS & PUNCHES



- Automatic self-holding
- Black color coded



351-100-550 Self-holding handle for forceps and punches (recommended handle)



Grasping Forceps with toothed cup, double action, 600mm working length

Insert only

Complete instrument

			Insert only	Complete instrument
 SEMI RIGID	semi-rigid	3 Charr., Ø 1.0mm	253-950-503	253-950-003
		4 Charr., Ø 1.3mm	253-950-504	253-950-004
		5 Charr., Ø 1.6mm	253-950-505	253-950-005
 FLEXIBLE	flexible	3 Charr., Ø 1.0mm	253-951-503	253-951-003
		4 Charr., Ø 1.3mm	253-951-504	253-951-004
		5 Charr., Ø 1.6mm	253-951-505	253-951-005



Stone Grasping Forceps for large stones,
double action jaws, 600mm working length

Insert only

Complete instrument



semi-rigid

3 Charr., Ø 1.0mm

253-952-503

253-952-003

4 Charr., Ø 1.3mm

253-952-504

253-952-004

5 Charr., Ø 1.6mm

253-952-505

253-952-005



flexible

3 Charr., Ø 1.0mm

253-953-503

253-953-003

4 Charr., Ø 1.3mm


253-953-504

253-953-004

5 Charr., Ø 1.6mm

253-953-505

253-953-005



HAMADA long Stone Grasping Forceps,
double action jaws, 600mm working length

Insert only

Complete instrument



semi-rigid

3 Charr., Ø 1.0mm

253-954-503

253-954-003

4 Charr., Ø 1.3mm

253-954-504

253-954-004

5 Charr., Ø 1.6mm

253-954-505

253-954-005



flexible

3 Charr., Ø 1.0mm

253-955-503

253-955-003

4 Charr., Ø 1.3mm

253-955-504

253-955-004

5 Charr., Ø 1.6mm

253-955-505

253-955-005



Biopsy Forceps, 600mm working length

Insert only

Complete instrument



semi-rigid

3 Charr., Ø 1.0mm

253-956-503

253-956-003

4 Charr., Ø 1.3mm

253-956-504

253-956-004

5 Charr., Ø 1.6mm

253-956-505

253-956-005



flexible

3 Charr., Ø 1.0mm

253-957-503

253-957-003

4 Charr., Ø 1.3mm

253-957-504

253-957-004

5 Charr., Ø 1.6mm

253-957-505

253-957-005


SCISSORS



- Standard handle
- White color coded




351-100-500 Standard handle for scissors
(recommended handle)



Scissors, single action, 600mm working length		Insert only	Complete instrument
 semi-rigid	3 Charr., Ø 1.0mm	253-958-503	253-958-003
	4 Charr., Ø 1.3mm	253-958-504	253-958-004
	5 Charr., Ø 1.6mm	253-958-505	253-958-005
 flexible	3 Charr., Ø 1.0mm	253-959-503	253-959-003
	4 Charr., Ø 1.3mm	253-959-504	253-959-004
	5 Charr., Ø 1.6mm	253-959-505	253-959-005


ASSEMBLY



STEP 1

Open the handle.
The ball socket pushes out.

Keep the handle open.



STEP 2

1 Opened jaw


2 Close the jaw - the ball head pushes out. Keep the jaw closed.



STEP 3

Push back the locking nut.


With both opened handle and jaws, connect the ball head and ball socket.



STEP 4

Insert the ball head into the socket.

Keep the handle opened.



STEP 5

Close the handle. Shaft and insert will be retracted.
Keep the handle closed.

Screw the locking nut onto the handle.





DONE

Check the connection for being tight.

The knob rotates 360° and the handle as well as the jaw opens and closes.

OPTICAL FORCEPS

■ for use with Ø 4mm Cystoscopes



Scope is **below** forceps



Scope is **above** Forceps



253-420-100*
Spoon Biopsy Forceps with round cup,
double action



253-426-100*
Mini Biopsy forceps with oval cup,
double action
for use with 30° cystoscope



253-421-100
Grasping Forceps with serrated jaws,
double action



253-431-100*
Biopsy forceps with **large** oval cup,
double action
for use with 30° cystoscope



253-422-100
Biopsy Punch Forceps, thru-cutting,
double action




253-425-100
Biopsy forceps with mini oval cup,
double action
for use with 4mm, 70° cystoscope



253-423-100
Scissors, hook type blades,
double action




253-427-100*
Stent Grasping forceps
with serrated jaws, double action
for removal of stents




253-429-100* Adapter for use of optical forceps with resectoscope sheaths


FLEXIBLE INSTRUMENTS



253-940-005 3-Prongs Grasper, flexible, 600mm working length, 5 Charr.




253-585-005 Dormia Metal Stone Basket, flexible, 4 wires, 600mm working length, 5 Charr.



Dormia Nitinol Stone Basket, 600mm working length, 5 Charr.


- 253-581-045 with 4 wires, complete with handle
- 253-582-045 Spare Basket only, 4 wires
- 253-582-055 Spare Basket only, 5 wires
- 253-582-065 Spare Basket only, 6 wires



Polyp Snare with teflon tube, 500mm working length with HF connection, 7 Charr.

253-580-007 with handle, complete

253-581-007 Spare Snare only



253-584-005 Zeiss Loop Catheter, package with 5 pieces, sterile, single use, 5 Charr.



FLEXIBLE ELECTRODES AND HF CABLES


Ball Electrode, 470mm working length

 253-370-003 3 Charr.

 253-370-005* 5 Charr.

 253-370-007 7 Charr.

 253-370-010 10 Charr.



Needle Electrode, 470mm working length



253-371-003 3 Charr.

253-371-005 5 Charr.

253-371-007 7 Charr.

Loop Electrode 470mm working length



253-372-005 5 Charr.

253-372-007 7 Charr.

Ball Electrode, 980mm working length



253-370-103 3 Charr.



Monopolar High Frequency cables, 4,5mm length


253-610-001* with 5 mm plug for Erbe HF-unit

253-610-002* with 4 mm plug for Martin HF-unit

253-610-003* with 8 mm plug for Valleylab HF-unit

253-610-004* with 4 mm plug for older HF-units


700-501-453* SDS plug



CYSTOSCOPY ACCESSORIES



301-701-100* Grey rubber cap for use on all working channels



253-990-800 French Catheter Scale




Cleaning Brush, flexible
253-961-103 Ø 3mm, 100mm length,
for use in 2.0 - 2.6mm working channels

253-961-102 Ø 2.5mm, 120mm length,
for 5 Charr. working channels




300-215-000
Cleaning Brush, Ø 13mm, 60mm brush length, 230mm overall length


RESECTOSCOPY



RZ RESECTOSCOPY SETS

Transurethral resection of the prostate (TURP) is a common surgery for BPH (Benign Prostate Hyperplasia). During TURP, the prostate is visualized through the urethra and tissue is removed using electric current. A monopolar or bipolar resectoscope is inserted transurethral and the enlarged prostatic tissue is removed with cutting and sealing electrodes.








Complete Instrument
consisting of **1** - **4**

Instrument components

- 1** Endoscope
- 2** Sheath system:
 - Outer Sheath
 - Inner Sheath
 - Obturator
- 3** Working Element
(active or passive)
- 4** Electrode



RZ **HYBRID** RESECTOSCOPY SYSTEM




Complete Instrument
consisting of **1 - 5**
- only cable needs to be changed


Instrument components

- 1** Endoscope
- 2** Hybrid Electrode
- 3** Sheath system:
 - Outer Sheath
 - Inner Sheath
 - Obturator
- 4** Working Element
(active or passive)
- 5** Monopolar / Bipolar
Cable





Complete Instrument
consisting of **1-5**
- only cable needs to be changed



ONE set
TWO techniques



Endoscopes, suitable for 24 Charr. System



0° DIRECTION OF VIEW
300mm working length
356mm overall length
253-804-000*



12° DIRECTION OF VIEW
300mm working length
356mm overall length
253-804-012*




30° DIRECTION OF VIEW
300mm working length
356mm overall length
253-804-030*




Monopolar working elements suitable for 24 Charr. System

- Lightweight, titanium handle parts
- Ergonomic handle design with finger rest and rotating titanium thumb ring





253-000-200* **Passive** Titanium Working Element




253-000-201* **Active** Titanium Working Element



Working elements suitable for HYBRID 24 Charr. system with bipolar and monopolar cable connection



253-000-400★ HYBRID **Passive** Titanium Working Element



253-000-401★ HYBRID **Active** Titanium Working Element




HYBRID



Sheaths, to be used with 24 Charr. system, monopolar and HYBRID


- 253-000-256* Continuous Flow Resectoscope Sheath
Rotating inner sheath for 360° TURP/TURB.
Inner shaft with non reflecting, oblique ceramic tip.
RZ QuickLock mechanism, with Obturator



- 253-000-266 Inner Shaft, 24 Charr., with oblique ceramic tip for 253-000-256



- 253-000-356 Continuous Flow Resectoscope Sheath, fixed,
inner sheath with non reflecting, oblique ceramic tip.
RZ QuickLock mechanism, with obturator




- 253-000-366 Inner Shaft with oblique ceramic tip, for 253-000-356



- 253-000-724 Central Valve Resectoscope Sheath, 24 Charr.,
oblique ceramic tip, with obturator



- 253-000-524 Resectoscope Sheath, 24 Charr.
oblique ceramic tip, with obturator



- 253-428-500 Adapter for resectoscope outer sheaths
253-000-256 and -356
and Ellik or bladder syringes



Obturator, to be used with 24 Charr. system, monopolar and HYBRID



253-000-044★ Standard Obturator, for 24/26 Charr. resectoscope sheaths



253-000-054★ Deflecting Obturator, for 24/26 Charr. resectoscope sheaths



253-000-084★ Visual Obturator for 4mm RZ Cystoscopes, for 24/26 Charr. resectoscope sheaths



253-000-094★ Visual Obturator with 5 Charr. working channel, for 4mm RZ Cystoscopes, for 24/26 Charr. resectoscope sheaths



Electrodes, to be used with
24 Charr. system



Cutting Loop, angled, for 12°/30° telescope

253-300-201 *



Cutting Loop, 90° angled, for 0° telescope

253-300-222



Cutting Loop, straight in longitudinal direction

253-300-202



Coagulation Electrode, pointed

253-300-203 *



Coagulation Electrode, conical

253-300-204



Coagulation Electrode, ball shape, Ø 3mm

253-300-205 *



Coagulation Electrode, ball shape, Ø 5mm

253-300-206 *



Coagulation Electrode, roller shape, Ø 3mm

253-300-207



Coagulation Electrode, roller shape, Ø 5mm

253-300-208



Electrodes, to be used with
24 Charr. system



Spike Electrode, Ø 3mm

253-300-210



Spike Electrode, Ø 5mm

253-300-209



Roller-Electrode, Ø 3mm

253-300-212



Roller-Electrode, Ø 5mm

253-300-211




Band Vapor Electrode, 1.2 mm width

253-300-213



Band Vapor Electrode, 90° angled

253-300-214



Band Vapor Electrode

253-300-215



Curette, sharp, for use without HF



253-300-218



Needle Electrode

253-300-219

24
CHARR.



Sterile disposable Electrodes,
to be used with 24 Charr. system



STERILE EO

Sales unit: package with 6 pieces,
each individually sterile packed.



Cutting Loop, angled, for 12°/30° telescope

253-300-201S*



Collins knife electrode

253-300-203S*




Ball electrode, Ø 3mm

253-300-205S*



Ball electrode, Ø 5mm

253-300-206S*



24
CHARR.

HYBRID



HYBRID Electrodes, to be used with 24 Charr. system

- Works in bipolar and monopolar mode
- Made of a standard wire of $\varnothing$ 0.35mm



Cutting Loop, angled, for 12°/30° telescope

253-400-200



Cutting Loop, angled, for 12°/30° telescope

253-400-201*



Cutting Loop, thick wire, angled, for 12°/30° telescope

253-400-202



Collins Coagulation Electrode, pointed

253-400-203*



Coagulation Electrode, conical

253-400-204



Coagulation Electrode, ball shape, $\varnothing$ 3mm

253-400-205*



Roller-Electrode, $\varnothing$ 3mm

253-400-211



Cutting Loop, 90° angled, for 0° telescope

253-400-222



Coagulation Electrode, ball shape, $\varnothing$ 5mm


253-400-206

24
CHARR.

HYBRID
HYBRID




HYBRID „Mushroom“ Electrodes, to be used with 24 Charr. system, only for bipolar use with NaCl




„Mushroom“ Electrode, blunt/flat

253-400-701




„Mushroom“ Electrode, blunt with slot, concave

253-400-702



„Mushroom“ Electrode,
half-round (convex), Ø 3mm, length 1.9mm

253-400-703



„Mushroom“ Electrode,
half-round (convex), Ø 3mm, length 2.5mm

253-400-704





XL RESECTOSCOPE


◀◀ 245mm working length ▶▶




253-804-530 XL Cystoscope, Ø 4mm, 30° direction of view, 365mm length



253-000-251 XL **active** working element, Titanium,
for use with 253-804-530 and 253-000-666



253-000-250 XL **passive** working element, Titanium,
for use with 253-804-530 and 253-000-666





253-000-666 XL Resectoscope Sheath, continuous flow, rotating,
for use with 253-000-251/-250 and 253-804-530

253-000-676 Outer Shaft






253-000-677 Inner Shaft


253-000-678 Obturator






XL 24 Charr. Monopolar Electrodes

	Cutting Loop, straight		253-300-252
	Cutting Loop, 15° angled		253-300-253
	Cutting Loop, 45° angled		253-300-255
	Cutting Loop, 60° angled		253-300-256
	Cutting Loop, 120° angled		253-300-251
	Roller Electrode		253-300-258
	Kife Electrode		253-300-259



Endoscopes, suitable for HYBRID 22 Charr. System



Ø
2.9 mm



● 0° DIRECTION OF VIEW

300mm working length
356mm overall length 351-829-000



○ 12° DIRECTION OF VIEW

300mm working length
356mm overall length 351-829-012





● 30° DIRECTION OF VIEW

300mm working length
356mm overall length 351-829-030




Monopolar Working Elements







Hybrid Working Elements and Sheath
suitable for HYBRID 22 Charr. System



253-000-410 **Passive** Titanium Working Element



253-000-411 **Active** Titanium Working Element



253-000-352 Resectoscope Sheath, rotating inner tube
for continuous suction / irrigation

253-000-555 Obturator only




Monopolar Electrodes, suitable for HYBRID 22 Charr. System

	Cutting Loop, angled	253-000-501
	Cutting Loop, straight	253-000-502
	Coagulation Electrode, Knife	253-000-503
	Coagulation Electrode, Ball	253-000-505
	Coagulation Electrode, Roller	253-000-507


Electrodes, suitable for HYBRID 22 Charr. System

		
	Cutting Loop, 30° angled	253-000-541
	Cutting loop, straight	253-000-542
	Knife Electrode	253-000-543
	Ball Electrode	253-000-545
	Roller Electrode, Ø 3mm	253-000-547

18.5
CHARR.



Endoscopes, suitable for HYBRID 18.5 Charr. System



Ø
2.9 mm



0° DIRECTION OF VIEW
300mm working length
356mm overall length 351-829-000



12° DIRECTION OF VIEW
300mm working length
356mm overall length 351-829-012




30° DIRECTION OF VIEW
300mm working length
356mm overall length 351-829-030




Working Elements and Sheaths,
suitable for HYBRID 18.5 Charr. System



351-400-200 **Passive Titanium Working Element**



351-400-201 **Active Titanium Working Element**



- 351-000-185 Resectoscope Continuous Flow Sheath,
with rotating sheath, incl. Obturator, with QuickLocking
- 351-000-155 Outer Shaft
- 351-000-195 Inner Shaft
- 351-000-186 Obturator



351-000-180 **Visual Obturator, with 5 Charr. working channel**



Electrodes, 16 Charr., suitable for HYBRID 18.5 Charr. System



Cutting Loop

351-000-201



Cutting loop, straight

351-000-202



Cutting loop, straight

351-000-203



Ball Electrode

351-000-205



Roller Electrode, 2mm

351-000-206



HIGH FREQUENCY CABLES

for RZ resectoscope systems, 4.5m length





253-610-001* with 5 mm plug for Erbe HF-unit

253-610-002* with 4 mm plug for Martin HF-unit

253-610-003* with 8 mm plug for Valleylab HF-unit

253-610-004* with 4 mm plug for older HF-units

700-501-453* SDS plug




253-400-601 for Erbe HF-unit


253-400-602 for Martin HF-unit

700-500-453 SDS plug

253-400-606 for BOWA



SUPRAPUBIC TROCAR SET



253-400-501 Trocar Sheath, 16 Charr., with pressure channel and spare silicone sealings



253-400-503 Suprapubic Trocar with pyramidal tip



253-400-504 Suprapubic Trocar with cannula tip




253-400-505 Aspiration Cannula



253-400-506 Overflow Cannula




253-400-508 Trocar guide sheath, lateral open



253-400-502 Suprapubic Trocar Stabilizing Plate, lateral open


LASER ENUCLEATION



LASER ENUCLEATION


Holmium laser prostate surgery is a minimally invasive treatment for Benign Prostate Hyperplasia (BPH). Laser enucleation of the prostate uses a laser resectoscope to enucleate the prostate tissue that is blocking urine flow. A morcellator which passes through the morcroscope is then used to cut the prostate and remove the adenoma.

As with other types of prostate laser surgery, Laser enucleation can offer faster recovery, less risk of bleeding and transfusion and better symptom relief compared to traditional prostate surgery.




Sheaths for Laser applications
from 26 to 18.5 Charr.

Front Fire Laser resection



LASER ENUCLEATION SETS FRONT FIRE

HoLEP / ThuLEP / Overview



24
CHARR.

22
CHARR.

18.5
CHARR.




253-804-000 0°
253-804-012 12°
253-804-030 30°




351-829-000 0°
351-829-012 12°
351-829-030 30°




351-829-000 0°
351-829-012 12°
351-829-030 30°




253-000-256



253-000-352



351-000-185



253-901-220
(to be used for morcellation)



253-000-300 fibers up to 0.8mm
253-000-312 fibers up to 1.2mm



253-000-319 fibers up to 0.8mm
253-000-316 fibers up to 1.2mm



253-000-318 fibers up to 0.8mm



253-000-305



253-900-121




253-905-220



253-905-220

24
CHARR.


HoLEP Set



Cystoscope, autoclavable, Ø 4.0mm

253-804-030 30° direction of view


253-804-012 12° direction of view



253-000-256* Continuous Flow 26 Charr. Resectoscope Sheath with rotating inner sheath and non-reflecting ceramic tip, including atraumatic obturator



253-000-266 Inner shaft, 24Charr., with oblique ceramic tip for 253-000-256



253-874-000 Laser Working Insert for use with Resectoscope Sheath and 4mm Cystoscope, with 7 Charr. working channel



RZ Laser Working Element, passive, for laser-probes

253-000-300 for fibres up to 0.8mm

253-000-312* for fibres up to 1.2mm



300-011-184* Luer-Lock Tuohy Borst Adapter (Ø 0.6 - 1.4mm)




RZ LL-Connection, rotatable for use with Laser Probe

253-000-301* short

253-000-302 long

24
CHARR.


HoLEP Set



Cystoscope, Ø 4.0mm

253-804-030 30° direction of view

253-804-012 12° direction of view



253-000-256* Continuous Flow 26 Charr. Resectoscope Sheath with rotating inner sheath and non-reflecting ceramic tip, including atraumatic obturator

253-000-305 RZ Laser Working Element, passive, with changeable guide for laser probes



1 Guide for Laser Probes

253-000-320 0.8mm

253-000-321 1.0mm

253-000-322 1.2mm

253-000-323 1.4mm

2 Guide for Laser Probes, with retractor

253-000-330 0.8mm


253-000-331 1.0mm

253-000-332 1.2mm

253-000-333 1.4mm



300-011-184* Luer-Lock Tuohy Borst Adapter (Ø 0.6 - 1.4mm)




(253-000-302)

RZ LL-Connection, rotatable
for use with Laser Probe

253-000-301* short

253-000-302 long




detachable


24
CHARR.

XL Laser Resection


◀ 245mm working length ▶



253-804-530 XL Cystoscope, Ø 4mm, 30° direction of view, 365mm length



253-000-666 XL Resectoscope Sheath, continuous flow, rotating, for use with 253-000-251/-250 and 253-804-530




253-000-350 XL Laser Working Element for fibers Ø1.2mm, passive



slim HoLEP[®]
Figueiredo


MiLEP Set acc. Dr. Felipe C. A. de Figueiredo



Cystoscope, Ø 2.9mm

253-829-030 30° direction of view

253-829-012 12° direction of view



253-000-352 Resectoscope Sheath 22 Charr., rotating inner tube
for continuous suction / irrigation

253-000-362 Inner tube only



253-000-082 Visual Obturator, 22 Charr.




RZ Laser Working Element, passive, for laser-probes

253-000-319 for fibres up to 0.8mm

253-000-316 for fibres up to 1.2mm



300-011-184[★] Luer-Lock Tuohy Borst Adapter (Ø 0.6 - 1.4mm)



(253-000-302)

RZ LL-Connection, rotatable
for use with Laser Probe


253-000-301[★] short

253-000-302 long


18.5
CHARR.

ultra slim HoLEP
figueiredo


MiLEP Set acc. Dr. Felipe C. A. de Figueiredo



Cystoscope, autoclavable, Ø 2.9mm
253-829-030 30° direction of view
253-829-012 12° direction of view




351-000-185 Resectoscope Continuous Flow Sheath, 18.5 Charr.,
with rotating sheath, incl. obturator, with QuickLocking
351-000-195 Inner tube only



253-000-318 RZ Laser Working Element, passive,
for laser-probes with fibres up to 0.8mm




300-011-184* Luer-Lock Tuohy Borst Adapter (Ø 0.6 - 1.4mm)




RZ LL-Connection, rotatable
for use with Laser Probe
253-000-301* short
253-000-302 long

24
CHARR.

TISSUE MORCELLATION



253-900-121* Nephroscope 6°, oval instrument channel, parallel eye piece, 250mm working length



253-000-256* Continuous Flow 26 Charr. Resectoscope Sheath with rotating inner sheath and non-reflecting ceramic tip, including atraumatic obturator




253-428-100* Adapter for use with RZ Nephroscopes and Resectoscope Outer sheaths 253-000-356/-256

22
CHARR.

slim HoLEP_κ Morcroscope

acc. Dr. Felipe C. A. de Figueiredo

figueiredo




253-905-220 Morcroscope, 22 Charr., 220mm working length,
5mm working channel*

18.5
CHARR.

* to be used for morcellation after laser enucleation with 18.5 Charr. Resectoscope

ultra slim HoLEP_κ



253-901-220 Continuous Flow Sheath, 22 Charr.,
220mm working length, 2 way stopcock,
quicklock, with Obturator




253-904-220 Sealing Cap for Morcroscope, with membrane to puncture through,
package of 10 pieces

Suitable sterilization trays can be found on page 120



UNIVERSAL LASER CYSTOSCOPE SET


Side fire, compatible with all laser fibres up to 2.5mm
for PVP, HoLEP & ThuVAP



Cystoscope, autoclavable, Ø 4.0mm

253-804-030 30° direction of view


253-804-012 12° direction of view



253-870-023★ High Continuous Flow Laser Cystoscope Sheath, 23 Charr.,
including sheath and blunt obturator




253-872-000★ Telescope bridge and laser guide with 1 working channel,
for 8 Charr. laser fibers and instruments




253-873-000 Visual Obturator for easy insertion of the sheath
into the urethra and bladder



LASER RESECTION SET for PVP, HoLEP & ThuVAP



Cystoscope, autoclavable, Ø 4.0mm
253-804-030 30° direction of view
253-804-012 12° direction of view




253-000-256* Continuous Flow 26 Charr. Resectoscope Sheath with rotating inner sheath and non-reflecting ceramic tip, including atraumatic obturator




253-874-000 Laser Working Insert with 7 Charr. working channel, for 24/26 Charr. sheaths and 4mm telescope

URETHROTOME




21
CHARR.

URETHROTOME SET




253-804-000* Cystoscope, 0°, Ø 4mm


253-000-200* Passive Titanium working element



253-500-221* Urethrotome Sheath oval shape, 21 Charr. with 5 Charr. working channel for catheters, bougies and accessories




253-500-321* Hollow Obturator with channel for guide wires



253-500-421 Visual Obturator, 21 Charr., for Ø 4mm telescope

253-500-230* Supplementary Sheath, half round, to insert Foley catheters, double Lock for secure fixation on the Urethrotome Sheath

253-500-325 Bridge with channel for 5 Charr. accessories



253-500-229 Continuous Flow supplementary sheath with double lock



1 Guide Stem Urethrotome Knives, for use without HF



253-500-561★

Sachse knife, straight



253-500-562

Sickle knife, retrograde



253-500-563★

Half round knife, 180° cutting




253-500-564

WaveCut knife, straight




253-300-999★ Storage container for knives and electrodes, metal, perforated



351-120-105 Flexible Biopsy forceps, 5 Charr., 40cm working length

15.5
CHARR.

URETHROTOME SET




253-829-000 RZ Cystoscope, Ø 2.9mm, 0° direction of view




253-500-500
Urethrotome Working Element
for use with RZ Cystoscope, Ø 2.9mm, 0°



253-500-125 Urethrotome Bridge 15.5 Charr.,
with 5 Charr. working channel,
for Ø 2.9mm endoscopes



253-500-115 15.5 Charr. Urethrotome Sheath with lateral channel
for guide wires and other accessories



253-500-315 Obturator for use with 253-500-115



253-500-581

Sachse knife, straight, 15 Charr.




253-500-582

Sickle knife, retrograde, 15 Charr.



253-500-100 OTIS Dilation Urethrotome complete with 3 atraumatic tips and 2 knives, in wooden case

253-500-103 Wooden case only



253-500-101 Urethrotome only





253-500-102 Knife only




253-500-108 Set with 3 atraumatic tips

LITHOTRIPSY






253-804-000* Cystoscope, 0°, Ø 4mm, autoclavable




253-755-020* Punch Sheath with Central Valve, 25 Charr. including Obturator




253-755-030 Visual Obturator with 7 Charr. working channel





253-755-040 Obturator only



253-755-010* Mauermayer Stone Punch




253-755-500* Hamada Bladder Stone Lithotripter, dismantable
for use with RZ Cystoscope 253-804-030 or 253-804-070




Comparison:
field of view

- 1 with 30° Endoscope
- 2 with 70° Endoscope




Dismountable




253-804-030 Cystoscope, autoclavable, Ø 4.0mm, 30° direction of view

To be used with Sheaths:




Cystoscope Sheath
253-300-025* 25 Charr. ○




Resectoscope Sheath
253-000-256* 26 Charr. ●
253-000-258 28 Charr. ●




Resectoscope Sheath
253-000-356 Fixed, 26 Charr. ●
253-000-358 Fixed, 28 Charr. ●



253-430-100* Adapter for use with Outer Resectoscope Sheaths




253-760-000* Stone Crushing Forceps for small stones, single action




253-804-030 Cystoscope, autoclavable, Ø 4.0mm, 30° direction of view


Suitable Resectoscope Sheaths



253-000-256* 26 Charr. ●
253-000-258 28 Charr. ●




253-000-356 Fixed, 26 Charr. ●
253-000-358 Fixed, 28 Charr. ●




253-780-000 Tissue Punch to cut up denatured tissue after Prostata Resection

NEPHROSCOPY





12 Charr. MINI NEPHROSCOPE



253-900-140* Mini Nephroscope, 12° direction of view, 220mm working length

12
CHARR.

Ø SHAFT

7.0
CHARR.

WORKING
CHANNEL*

* for use with instruments up to 5 Charr.



253-905-400
Silicone Valve,
package with 5 pieces



253-900-010
Silicone Valve, package with 10 pieces,
for use with very delicate instruments
from 0.5 – 1.6mm diameter.



OPERATING SHEATHS for use with RZ Mini Nephroscope


Mini PCNL Sheaths with valve, 130mm working length



253-900-141★ 15.0 Charr. ●




253-900-142★ 16.5 Charr. ●



253-900-143★ 21.0 Charr. ●

Long Mini PCNL Sheaths, 180mm working length



253-900-171 15.0 Charr.



253-900-175 16.5 Charr.



253-900-173 21.0 Charr.

One-Step Dilators with 1 central guide wire channel




253-900-146★ for 253-900-141 + 253-900-171



253-900-147★ for 253-900-142 + 253-900-175, with 1 excentric guide wire channel




253-900-148★ for 253-900-143 + 253-900-173, with 1 excentric guide wire channel




253-900-145★ 2-piece Sealing Adapter with sheath
and obturator for use with 253-900-141

RZ NEPHROSCOPY INSTRUMENTS


semi-rigid, 5 Charr.



Stone Grasping Forceps,
double action jaws, retrograde serration
253-950-455* 450mm working length




Stone Grasping Forceps
for large stones, double action jaws, 1 x 1 teeth
253-952-455* 450mm working length




Hamada Long Stone Grasping Forceps,
double action jaws, fine serrated
253-954-455 450mm working length




Biopsy Forceps,
double action jaws
253-956-455* 450mm working length



Scissors,
single action jaws
253-958-455 450mm working length



253-960-360* Stone Grasping Forceps, rigid, 5 Charr., 400mm,
for grasping of bigger stones and -fragments
triple toothed jaw, U-spring handle




253-900-200 Sterile disposable Guide wire, J-tip, flexible, 150cm length



RZ COMPACT NEPHROSCOPE

220mm working length




Compact Nephroscope, including Irrigation-Valve-Adapter 253-900-000

18.8 Charr., 12.9 Charr working channel

253-900-815 15°

253-900-830 30°



Suitable Shaft, with obturator, 23 Charr., 220mm working length


253-900-888



253-900-000

Irrigation-Valve-Adapter

Percutane Nephroscopy, 190mm length



RZ Nephroscope, oblique view, with 2 Luer-Lock stopcocks and detachable instrument stopcock

253-900-100 6° direction of view

253-900-150 15° direction of view




253-900-105 Stopcock, only for use with RZ Nephroscopes



RZ Nephroscope, parallel view, with 2 Luer-Lock stopcocks and detachable instrument stopcock

253-900-101 6° direction of view

253-900-151 15° direction of view



Operating Sheath

1 253-900-110* 26 Charr.

2 253-900-111 24 Charr.



Hollow Obturator

1 253-900-115* 26 Charr.

2 253-900-116 24 Charr.



Percutane Nephroscopy, 250mm length



RZ Nephroscope, oblique view, with 2 Luer-Lock stopcocks and detachable instrument stopcock, 250mm working length

253-900-120 6° direction of view

253-900-160 15° direction of view



253-900-105 Stopcock, only for use with RZ Nephroscopes (including: 253-900-100 and 253-900-101)



Nephroscope Scope, oval instrument channel, parallel eye piece, 250mm working length

253-900-121 6° direction of view

253-900-161 15° direction of view



Operating Sheath

1 253-900-130* 26 Charr.


2 253-900-131 24 Charr.



Hollow Obturator

1 253-900-135* 26 Charr.

2 253-900-136 24 Charr.






253-900-155★ Amplatz Metal Dilator set containing 6 Dilators (9, 12, 15, 18, 21, 24 Charr.)
2 flexible and 2 rigid guide rods for use with guide wires



Amplatz Metal Dilator
253-900-156★ 27 Charr.
253-900-157★ 30 Charr.



253-900-201★ Puncture Cannula with inner and outer cannula




253-900-202★ Dilation Cannula, 3mm Ø with inner and outer cannula



253-900-200 Sterile disposable Guide wire, J-tip, flexible, 150cm length





Stone Grasping forceps with fenestrated jaws,
double action, with flushing port
253-900-210* 360mm working length
253-900-250* 430mm working length



Stone Grasping forceps with coarse serrated jaws,
double action
253-900-212* 360mm working length
253-900-252* 430mm working length




Biopsy forceps with sharp cup,
single action
253-900-211 360mm working length
253-900-251 430mm working length



Scissors, sharp tip,
single action
253-900-213 360mm working length
253-900-253 430mm working length

Stone Grasping Forceps with spring handle and finger platform
Ø 4mm, Luer-Lock connection




Fenestrated jaws
253-900-205 360mm
253-900-245 450mm



3 prongs
253-900-206 360mm
253-900-246 450mm

URETERO-RENOSCOPY



RZ URETERO-RENOSCOPES


- FullHD 50.000 pixel fiber system
- Semi-flexible shaft
- Silicone valve system with 2 parallel sealings
- Large working channels
- Irrigation/suction channels



Sealing caps included
301-701-060 Ø 0.6mm
301-701-030 Ø 0.3mm



Atraumatic tip design



Tapered shaft



90° stopcocks



5° direction of view
65° field of view

ADVANCED SEMI-RIGID URETERO-RENOSCOPES

standard LINE

253-905-540* 425mm
253-905-530 310mm

7.5
CHARR.

DISTAL TIP

9.0
CHARR.

DISTAL

13
CHARR.

PROXIMAL

5.8
CHARR.

WORKING
CHANNEL*

* for use with instruments up to 5 Charr.

slim LINE

253-910-540* 425mm
253-910-530 310mm

6.0
CHARR.

DISTAL TIP

7.5
CHARR.

DISTAL

10.5
CHARR.

PROXIMAL

5.0
CHARR.

WORKING
CHANNEL*

* for use with instruments up to 4 Charr.

superslim LINE

253-915-540* 425mm
253-915-530 310mm

5.0
CHARR.

DISTAL TIP

6.5
CHARR.

DISTAL

8.0
CHARR.


PROXIMAL

4.0
CHARR.


WORKING
CHANNEL*

* for use with instruments up to 3 Charr.


VERSATIP 360° rotating, color coded Instruments




351-100-550
Self-holding handle for
forceps and punches
(recommended handle)



351-100-500
Standard handle
for scissors
(recommended handle)



Grasping Forceps with toothed cup,
double action, 600mm working length




		Insert only	Complete instrument
 SEMI RIGID	semi-rigid	3 Charr., Ø 1.0mm	253-950-503
		4 Charr., Ø 1.3mm	253-950-504
		5 Charr., Ø 1.6mm	253-950-505
 FLEXIBLE	flexible	3 Charr., Ø 1.0mm	253-951-503
		4 Charr., Ø 1.3mm	253-951-504
		5 Charr., Ø 1.6mm	253-951-505



Stone Grasping Forceps for large stones,
double action jaws, 600mm working length




		Insert only	Complete instrument
 SEMI RIGID	semi-rigid	3 Charr., Ø 1.0mm	253-952-503
		4 Charr., Ø 1.3mm	253-952-504
		5 Charr., Ø 1.6mm	253-952-505
 FLEXIBLE	flexible	3 Charr., Ø 1.0mm	253-953-503
		4 Charr., Ø 1.3mm	253-953-504
		5 Charr., Ø 1.6mm	253-953-505



HAMADA long Stone Grasping Forceps, double action jaws, 600mm working length		Insert only	Complete instrument
 semi-rigid	3 Charr., Ø 1.0mm	253-954-503	253-954-003
	4 Charr., Ø 1.3mm	253-954-504	253-954-004
	5 Charr., Ø 1.6mm	253-954-505	253-954-005
 flexible	3 Charr., Ø 1.0mm	253-955-503	253-955-003
	4 Charr., Ø 1.3mm	253-955-504	253-955-004
	5 Charr., Ø 1.6mm	253-955-505	253-955-005




Biopsy Forceps, 600mm working length		Insert only	Complete instrument
 semi-rigid	3 Charr., Ø 1.0mm	253-956-503	253-956-003
	4 Charr., Ø 1.3mm	253-956-504	253-956-004
	5 Charr., Ø 1.6mm	253-956-505	253-956-005
 flexible	3 Charr., Ø 1.0mm	253-957-503	253-957-003
	4 Charr., Ø 1.3mm	253-957-504	253-957-004
	5 Charr., Ø 1.6mm	253-957-505	253-957-005



Scissors, single action, 600mm working length		Insert only	Complete instrument
 semi-rigid	3 Charr., Ø 1.0mm	253-958-503	253-958-003
	4 Charr., Ø 1.3mm	253-958-504	253-958-004
	5 Charr., Ø 1.6mm	253-958-505	253-958-005
 flexible	3 Charr., Ø 1.0mm	253-959-503	253-959-003
	4 Charr., Ø 1.3mm	253-959-504	253-959-004
	5 Charr., Ø 1.6mm	253-959-505	253-959-005

FLEXIBLE INSTRUMENTS



253-960-025 Stone Basket, 2.5 Charr., sterile, single use, package of 5 pieces



253-962-003 Brush for Cytology, 3.0 Charr., for single use, package of 5 pieces



253-940-005 3-Prongs Grasper, flexible, 600mm working length, 5 Charr.



253-585-005 Dormia Metal Stone Basket, flexible, 4 wires, 600mm working length, 5 Charr.


Accessories




253-905-400
Silicone Valve,
package with 5 pieces




253-900-010
Silicone Valve, package with 10 pieces,
for use with very delicate instruments
from 0.5 – 1.6mm diameter.



600-662-077
Endoscopic wire basket for URS
640 x 150 x 77mm




253-905-000
Container for URS with silicone inserts
640 x 140 x 60mm




Tray for URS and instruments
600-662-540 Tray
600-600-191 XL Container, Standard,
bottom non perforated,
732 x 280 x 190mm,
silver lid (*)


PEDIATRIC




CYSTO-URETHROSCOPE SET FOR NEW BORN



RZ Cystoscope, Ø 2.0mm
253-221-030 30°
253-221-000 0°





253-200-007 Cysto-Urethroscope Sheath, 7.5 Charr.



253-200-009 Cysto-Urethroscope Sheath, 9 Charr.,
with 3 Charr. Instrument Channel

3 Charr. Instruments for new born set




253-202-003
Alligator Grasping Forceps
280mm working length



253-204-003
Biopsy Forceps,
280mm working length




253-370-253
Ball Electrode,
250mm working length




253-371-253
Needle Electrode,
250mm working length


Ø 2.0mm PEDIATRIC CYSTOSCOPY




RZ Pediatric Fiber Cystoscope Ø 2.0mm, 187.5mm length
 253-219-100 0°
 253-219-130 30°
 253-219-170 70°



253-200-075 Pediatric Cystoscope Sheath, 7.5 Charr.




253-200-095 Pediatric Cystoscope Sheath, 9.5 Charr.
 with 4 Charr. instrument channel



253-200-011 Pediatric Cystoscope Sheath, 11 Charr.
 with 5 Charr. instrument channel

Ø 2.7mm PEDIATRIC CYSTOSCOPY




RZ Pediatric Cystoscope, 2.7mm diameter, 187.5mm working length


253-227-000 0°

253-227-030 30°


253-227-070 70°



253-220-011 Pediatric Cystoscope Sheath, 11 Charr.
with 3 Charr. working channel




253-220-014 Pediatric Cystoscope Sheath 13 Charr.
with 3 Charr. working channel, for use
with Bridge or Albarran deflecting mechanism



253-220-025
Albarran deflecting mechanism for use with 253-220-014

253-220-030
Bridge for use with 253-220-014

3 Charr. Instruments for Ø 2.0 / 2.7mm Pediatric Cystoscopy




Ø 2.9mm ADOLESCENT CYSTOSCOPY



RZ Cystoscope, Ø 2.9mm
253-829-000 0°
253-829-030 30°



253-310-014 Cysto-Urethroscope Sheath 15 Charr., with Obturator and integrated working channel for 5 Charr. flexible instruments and 2 lateral Luer-Lock stopcocks ●



253-310-015 Cysto-Urethroscope Sheath 16 Charr., with Obturator and 2 lateral Luer-Lock stopcocks ●




1 253-310-100 Telescope bridge without channel
2 253-310-101 Telescope bridge with 1 channel for use with 253-310-014 and 5 Charr. flexible instruments




253-310-107 Albarran deflecting mechanism with 1 channel for 253-310-014, for use with one 5 Charr. instrument

Ø 2.0mm PEDIATRIC RESECTOSCOPE




253-250-200 Passive working element for use with Ø 2.0mm
RZ Fiber Cystoscope 253-219-100

RZ Pediatric Fiber Cystoscope Ø 2.0mm, 187.5mm length
253-219-100 0°
253-219-130 30°



253-250-211 Resectoscope Sheath, 11 Charr.
with obturator ●




253-250-213 Resectoscope Sheath, 13 Charr.
with obturator ●

Electrodes

11
CHARR.

13
CHARR.



Cutting Loop

253-251-111

253-252-111

Coagulation Electrode, blunt

253-251-112

253-252-112

Hook Electrode with ball

253-251-113

253-252-113

Hook Electrode without ball

253-251-114


253-252-114

Needle Electrode, pointed



253-251-115

253-252-115


Ø 2.0mm PEDIATRIC URETHROTOME



RZ Pediatric Fiber Cystoscope Ø 2.0mm, 187.5mm working length
 253-219-100 0°
 253-219-130 30°



253-250-200 Passive working element for use with Ø 2.0mm
 RZ Fiber Cystoscope 253-219-100



253-265-010
 Pediatric Urethrotome sheath
 10 Charr. with obturator



253-255-001 Bridge with 1 channel

Electrodes



Knife with straight cutting edge 253-265-101



Round knife 253-265-102



Sickle Knife 253-265-103



Hook Knife 253-265-104

RZ PEDIATRIC CYSTO-URETHRO-RENOSCOPE

standard LINE

253-200-100 150mm

7.5 CHARR.	9.5 CHARR.	6.0 CHARR.
DISTAL	PROXIMAL	WORKING CHANNEL*

* for use with instruments up to 5 Charr.

slim LINE

253-200-500 150mm

6.0 CHARR.	7.5 CHARR.	5.0 CHARR.
DISTAL	PROXIMAL	WORKING CHANNEL*




* for use with instruments up to 4 Charr.

superslim LINE


253-200-600 150mm

4.5 CHARR.	6.0 CHARR.	3.5 CHARR.
DISTAL	PROXIMAL	WORKING CHANNEL*


* for use with instruments up to 3 Charr.




INSTRUMENTS FOR PEDIATRIC CYSTO-URETHROTOMES



	253-200-128	253-204-128	253-203-128	Stone Grasping Forceps, double action jaws
	253-200-228	253-204-228	253-203-228	Large Stone Grasping Forceps, double action jaws
	253-200-328	253-204-328	253-203-328	Biopsy Forceps, double action
	253-200-428	253-204-428	253-203-428	Scissors, single action



	253-370-005	253-370-004	253-370-003	Ball Electrode
	253-371-005		253-371-003	Needle Electrode



HIGH FREQUENCY CABLES

for RZ pediatric resectoscope systems, 4.5m length





253-610-001★ with 5 mm plug for Erbe HF-unit

253-610-002★ with 4 mm plug for Martin HF-unit

253-610-003★ with 8 mm plug for Valleylab HF-unit

253-610-004★ with 4 mm plug for older HF-units

700-501-453★ SDS plug




253-400-601 for Erbe HF-unit


253-400-602 for Martin HF-unit


700-500-453 SDS plug

253-400-606 for BOWA



DEVICES AND FLEXIBLE INSTRUMENTS





253-800-000 Flexible Cystoscope


Shaft size: 15.5 Charr.
Working Channel: $\varnothing$ 7 Charr.
Angulation up/down: 210°/120°
Working length: 380mm



Shaft size: 15.5 Charr.
Slimline Sheath for
easy insertion



Angulation: 210°/120°



7 Charr. working
channel



Universal suction system



253-800-900 Flexible Uretero-Renoscope


Shaft size: 7.5 Charr.
WorkingChannel: $\varnothing$ 3.6 Charr.
Angulation up/down: 275°/275°
Working length: 680mm



Shaft size: 7.5 Charr.
Slimline sheath for
easy insertion



Angulation:
275° / 275°




3.6 Charr. working
channel



Ergonomic grip

ACCESSORIES



SYRINGES & EVACUATORS


with Turn Lock connection to connect to RZ Cystoscope and Resectoscope sheaths



TOOMEY Bladder Syringe
253-000-850 50cc



Bladder Syringe		Spare Glass Cylinder	
253-000-975	75cc	253-000-976	75cc
253-000-900	100cc	253-000-901	100cc
253-000-950*	150cc	253-000-951	150cc



Ellik Evacuator
253-500-000 with turn lock connection
253-500-001 with standard connection

- 1 253-500-003 Rubber ball
- 2 253-500-004 Glass body
- 3 253-500-020 Turn Lock connector with silicone tube
- 4 253-500-010 Standard connector with silicone tube
- 5 253-500-030 Quick Lock connection with plastic tube

ADAPTORS & VALVES




300-011-181* Luer-Lock tube connector,
9mm



300-011-506* Luer-Lock tube connector,
6mm



300-011-184* Luer-Lock Tuohy
Borst Adapter



300-011-178 Luer-Lock one-way stop-cock
with spring screw cap



300-011-179* Luer-Lock tube connector
with stopcock



300-011-515* Luer-Lock one-way
stopcock




Rubber sealing cap
301-701-100* Grey / Ø 1.0mm
301-701-060 Green / Ø 0.6mm
301-701-030 Red / Ø 0.3mm



351-905-105 Disposable Silicone Valve
with Luer-Lock connection
for use with 3 - 10 Charr.
instruments, package with
10 pieces, sterile packed



300-400-900* Rubber Seal for
Luer-Lock connectors



300-400-901* Rubber Seal for Luer-Lock
connectors

CATHETERS



Guyon Catheter Spanner & Introducer

- 1** 253-310-390* curved, 390mm
- 2** 253-310-392 straight, 400mm



320-603-455 Foley Catheter Spanner & Introducer, curved, 455mm




Catheter for male

320-720-080	8 Charr.	320-720-180	18 Charr.
320-720-100	10 Charr.	320-720-200	20 Charr.
320-720-120	12 Charr.	320-720-220	22 Charr.
320-720-140	14 Charr.	320-720-240	24 Charr.
320-720-160	16 Charr.		



Catheter for female

320-721-080	8 Charr.
320-721-100	10 Charr.
320-721-120	12 Charr.
320-721-140	14 Charr.
320-721-160	16 Charr.
320-721-180	18 Charr.
320-721-200	20 Charr.
320-721-220	22 Charr.
320-721-240	24 Charr.



Rotunda Catheter, s-curved


320-723-080	8 Charr.
320-723-100	10 Charr.
320-723-120	12 Charr.
320-723-140	14 Charr.
320-723-160	16 Charr.
320-723-180	18 Charr.
320-723-200	20 Charr.
320-723-220	22 Charr.
320-723-240	24 Charr.

BOUGIES



Guyon Bougie, strongly curved

253-340-335*	8 Charr.
253-341-335	9 Charr.
253-342-335*	10 Charr.
253-343-335	11 Charr.
253-344-335*	12 Charr.
253-345-335*	14 Charr.
253-346-335*	16 Charr.
253-347-335*	18 Charr.
253-348-335*	20 Charr.




253-349-335*	22 Charr.	253-354-335*	28 Charr.
253-350-335*	24 Charr.	253-355-335*	30 Charr.
253-351-335	25 Charr.	253-356-335	32 Charr.
253-352-335*	26 Charr.	253-357-335	34 Charr.
253-353-335	27 Charr.		



Van Buren Bougie, acutely curved

253-650-010*	10 Charr.
253-650-012*	12 Charr.
253-650-014*	14 Charr.
253-650-016*	16 Charr.
253-650-018*	18 Charr.
253-650-020*	20 Charr.
253-650-022*	22 Charr.


253-650-024*	24 Charr.
253-650-026*	26 Charr.
253-650-028*	28 Charr.
253-650-030*	30 Charr.
253-650-032	32 Charr.
253-650-034	34 Charr.
253-650-036	36 Charr.



Dittel Bougie, slightly curved

253-320-335*	8 Charr.
253-321-335	9 Charr.
253-322-335*	10 Charr.
253-323-335	11 Charr.
253-324-335*	12 Charr.
253-325-335*	14 Charr.
253-326-335*	16 Charr.
253-327-335*	18 Charr.
253-328-335*	20 Charr.

253-329-335*	22 Charr.
253-330-335*	24 Charr.
253-331-335	25 Charr.
253-332-335*	26 Charr.
253-333-335	27 Charr.
253-334-335*	28 Charr.
253-335-335*	30 Charr.
253-336-335	32 Charr.
253-337-335	34 Charr.



BOUGIES & DILATORS



OTIS Bougie à Boule

253-370-009 9 - 12 - 15 - 18 Charr.

253-370-021 21 - 24 - 27 - 30 Charr.



Dittel Bougie, straight

253-320-225	8 Charr.	253-329-225	22 Charr.
253-321-225	9 Charr.	253-330-225	24 Charr.
253-322-225	10 Charr.	253-331-225	25 Charr.
253-323-225	11 Charr.	253-332-225	26 Charr.
253-324-225	12 Charr.	253-333-225	27 Charr.
253-325-225	14 Charr.	253-334-225	28 Charr.
253-326-225	16 Charr.	253-335-225	30 Charr.
253-327-225	18 Charr.	253-336-225	32 Charr.
253-328-225	20 Charr.	253-337-225	34 Charr.



RZ Orifice Dilator

253-375-015 15 - 30 Charr.



Ureteral Dilator, flexible, 500mm

253-920-090 9.0 Charr.

253-920-135 13.5 Charr.

253-920-105 10.5 Charr.

253-920-150 15.0 Charr.

253-920-120 12.0 Charr.




320-998-000 Stockmann Meatus Clamp, 8cm



320-999-000* Strauss Penis Seizing Forceps, 13cm

CLEANING




SAFE + CLEAN ENDOSCOPE PLASTIC CONTAINER

with silicone bars on base and lid


600-505-000

Endoscopes	1 (Ø 1-10mm)
Endoscope length	max. 270mm
Dimensions	271 x 70 x 35mm




600-505-002

Endoscopes	max. 2
Endoscope length	max. 190mm
Dimensions	271 x 70 x 35mm



600-506-002*

Endoscopes	max. 2
Endoscope length	max. 390mm
Dimensions	394 x 70 x 35mm




600-501-000

Endoscopes 1

600-501-002


Endoscopes max. 2

Endoscope length	max. 450mm
Dimensions	457 x 70 x 35mm




ENDOSCOPE & INSTRUMENT CONTAINERS

600-508-000	for scopes and light guide cables
Endoscopes	max. 3 (Ø 5-10mm)
Fiber Optic Light Cable	1
Endoscope length	max. 450mm
Dimensions	465 x 212 x 51mm



600-464-229	for fixation of instruments with rubber ties
Dimensions	464 x 229 x 51mm




600-511-000	for instruments with rubber ties on 2 levels
Dimensions	502 x 239 x 102mm




ENDOSCOPE CLEANING BASKETS & CONTAINERS

with silicone bars on base and lid




600-292-052*

Endoscopes	max. 2
Endoscope length	max. 190mm
Dimensions	290 x 80 x 52mm



600-462-052*

Endoscopes	max. 2
Endoscope length	max. 350mm
Dimensions	460 x 80 x 52mm




600-460-055

Endoscopes 1

600-460-056


Endoscopes max. 2

Endoscope length	max. 340mm
Dimensions	455 x 86 x 43mm




600-450-070

Endoscopes	max. 2
Endoscope length	max. 330mm
Dimensions	450 x 70 x 70mm



600-662-077


Endoscopic wire basket for URS
640 x 150 x 77mm





253-905-000

Container for URS with silicone inserts
640 x 140 x 60mm


TRAYS FOR MiLEP SYSTEM



600-605-220 MiLEP Wire Basket
480 x 300 x 60mm



600-605-411 Wire basket for 22 Charr. Resectoscope
480 x 125 x 60mm





1/1


Length

Stainless Steel Basket,
round perforation, with lid

545mm / 485mm

complete

600-540-535	545 x 255 x 35mm
600-540-555	545 x 255 x 55mm
600-540-575	545 x 255 x 75mm
600-540-605	545 x 255 x 105mm
600-485-535	485 x 255 x 35mm
600-485-555	485 x 255 x 55mm
600-485-575	485 x 255 x 75mm
600-485-605	485 x 255 x 105mm




Lid only

600-540-500	545 x 255 mm
600-485-500	485 x 255 mm




Bottom only

600-540-530	545 x 255 x 30mm
600-540-550	545 x 255 x 50mm
600-540-570	545 x 255 x 70mm
600-540-600	545 x 255 x 100mm
600-485-530	485 x 255 x 30mm
600-485-550	485 x 255 x 50mm
600-485-570	485 x 255 x 70mm
600-485-600	485 x 255 x 100mm



Silicone Mat, green

600-281-004	520 x 230mm
600-281-003	440 x 230mm




3/4

Stainless Steel Basket,
round perforation, with lid


Length

405mm


600-405-535 405 x 255 x 35mm
600-405-555 405 x 255 x 55mm
600-405-575 405 x 255 x 75mm
600-405-605 405 x 255 x 105mm




Lid only
600-405-500 405 x 255 mm



Bottom only
600-405-530 405 x 255 x 30mm
600-405-550 405 x 255 x 50mm
600-405-570 405 x 255 x 70mm
600-405-600 405 x 255 x 100mm



600-281-002 Silicone Mat green, 380 x 230mm




1/2

Stainless Steel Basket,
round perforation, with lid

Length

255mm

600-255-535	255 x 245 x 35mm
600-255-555	255 x 245 x 55mm
600-255-575	255 x 245 x 75mm
600-255-605	255 x 245 x 105mm




Lid only

600-255-500 255 x 245 mm



Bottom only

600-255-530	255 x 245 x 30mm
600-255-550	255 x 245 x 50mm
600-255-565	255 x 245 x 70mm
600-255-600	255 x 245 x 100mm



600-281-001 Silicone Mat green, 220 x 230mm

Cleaning



600-100-000
RZ-Aquajet
Water Jet Pistol complete
with 9 tips and 1m pressure
tubing, with strap and
suspension hook




600-050-000
Hand brush, autoclavable
(10 pcs. per box)





600-950-100
RZ - Brush dispenser
for 10 brushes with suspension hook



700-900-220 Special Cotton for polishing metal sheaths and instruments



660-166-100
Holding Device for cleaning & sterilization



700-900-160
Instrument Oil in 50ml bottle



700-900-250
Oil Dispenser, 140mm



Cleaning brush (1 piece per box), 500mm length
300-135-905 for Ø 5 mm
300-135-910 for Ø 10 mm



300-215-000 Cleaning Brush, Ø 13mm, 215mm length, 60mm brush length



300-420-065 Cleaning Brush, Ø 3mm, 420mm length, 65mm brush length

EQUIPMENT



FLUID MANAGEMENT



Features:

- Versatile Irrigation Pump
- Pump leveling adjustment
- Pressure controlled roller pump

DELIVERY CONTAINS:

- Pump device
- Tubing set (253-910-001)
- Mains cable
- 2 Spare fuses
- Operating manual


253-910-000 Basic Irrigation Pump P1

Technical data

Line voltage	100-240 VAC
Power input	80VA
Irrigation power	0 - 1.0L/min
Pressure range	10 - 150mmHG
Measurements	380 x 130 x 400mm (WxHxD)
Weight	approx. 11kg

Accessories

253-910-001	Irrigation Tubing Set, silicone, autoclavable	non-sterile	included
253-910-002	Irrigation Tubing Set, PVC, 20 pieces, daily tube	sterile	optional
253-910-003	Patient tube, PVC, 2.5m, 40 pieces	sterile	optional



4K CAMERA SYSTEM



300-001-840 4K 3-Chip UHD-Camera K 4000

4K

Features:

- 4K UHD resolution with 3-chip technology
- Image algorithm enhancement of color and structure
- 4K video recording
- 4K image saving
- Camera head with 4 programmable buttons
- Boost function
- Fog reduction
- Freeze function
- 2.5x digital zoom
- 10 USER mode configurations
- Multilingual on screen menu
- SDI 4K Output
- Text generator (optional with keyboard)

DELIVERY CONTAINS:

- Main Unit (300-001-840)
- Camera head (300-001-841)
- HDMI cable
- USB Pen Drive
- Mains cable
- Operating manual

Technical data

Line voltage	100-240 VAC
Power consumption	60VA
Image sensor	3 MOS
Resolution	3840 x 2160 pixels (4K)
Index ratio	16:9
N/N ratio	54db

Video outputs	1x HDMI 4K 1x HDMI FullHD 1x SDI 4K 1x SDI FullHD
Dimensions	Control unit: 355 x 103 x 350mm (WxHxD) Weight control unit: 6.2kg Weight camera head: 145g

Accessories

300-001-841
UHD 3-Chip Camera Head
(included in the set)300-001-615
Title generator keyboard, waterproof IP68
US-Layout, cable 1.8m, color: grey,
305x140x15mm, 600gr.300-001-610
Footswitch for HD Cameras
cable length: 2m

FULL HD CAMERA SYSTEM



Features:

- Full HD camera with 3-chip MOS technology
- Pure Natural Colors technology
- Full HD video recording (F4 only)
- Full HD image saving (F4 only)
- Camera head with 3 programmable buttons
- Boost function
- Freeze function
- 2.5x digital zoom
- 3 USER mode configurations
- Multilingual on screen menu
- Different digital video outputs
- Text generator (optional with keyboard)

DELIVERY CONTAINS:

- Main Unit
- Camera head
- 2x BNC cable
- HDMI - DVI cable
- HDMI cable
- S-Video cable (F3 only)
- USB Pen Drive (F4 only)
- Mains cable
- Operating manual

300-001-800 3 Chip HD-Camera F3
300-001-804 3 Chip HD-Camera F4 (with recording function)

FULL HD

Technical data

Line voltage	100-240 VAC
Power input	20W
Image sensor	3 MOS
Resolution	1920 x 1080 pixels (FullHD)
Image resolution	>1000 TV lines, horizontal
Index ratio	16:9
Min. sensitivity	< 3lx (1/1.4)
Shutter speed	1/50 (1/60) to 1/100.000sec

Video outputs (F4)

2 x HD SDI
1x HDMI / DVI 2x CVBS
1 x YPbPr / RGB
1 x Y/C (Pal/NTSC)
1x CVBS (Pal/NTSC)

Video outputs (F3)

1 x HDMI / DVI
1 x YPbPr (Component)
1 x Y/C (Pal/NTSC)
1x CVBS (Pal/NTSC)

Dimensions

Control unit: 380 x 120 x 360mm (WxHxD)
Weight control unit: 11kg
Weight camera head: 145g

Accessories



300-001-802
High-End 3-MOS Camera Head
FullHD (included in the set)



300-001-615
Title generator keyboard, waterproof IP68
US-Layout, cable 1.8m, color: grey,
305x140x15mm, 600gr.



300-001-610
Footswitch for HD Cameras
cable length: 2m

FULL HD CAMERA SYSTEM



Features:

- Full HD camera with 1-chip MOS technology
- Natural color technology
- Full HD video recording (F2 Pro only)
- Full HD image saving (F2 Pro only)
- Camera head with 3 programmable buttons
- Boost function
- Freeze function
- 2.5x digital zoom
- 3 USER mode configurations
- Multilingual on screen menu
- HDMI video outputs
- Text generator (optional with keyboard)

300-001-600 HD-Camera F1
300-001-601 HD-Camera F2 Pro (with recording function)

FULL HD

Technical data

Line voltage	100-240 VAC
Power input	F2 Pro: 30W / F1: 20W
Image sensor	1/3" MOS
Resolution	1920 x 1080 pixels (FullHD)
Image resolution	>1000 TV lines, horizontal
Index ratio	16:9
Min. sensitivity	< 3lx (1/1.4)

DELIVERY CONTAINS:

- Main Unit
- Camera head
- 2x BNC cable
- HDMI - DVI cable
- HDMI cable
- S-Video cable
- USB Pen Drive (F2 Pro only)
- Operating manual
- Mains cable

Video outputs (F2 Pro)

2 x HD SDI (optional)
1 x HDMI / DVI
1 x Y-Pb-Pr (Component)
1 x Video (PAL/NTSC)
1x S-Video (PAL/NTSC)

Video outputs (F1)

1 x HDMI
1 x YPbPr / RGB
1 x S-Video (Y/C; Pal / NTSC)
1 x Video (CVBS; Pal / NTSC)

Dimensions

Control unit: 380 x 120 x 360mm (WxHxD)
Weight control unit: 11kg
Weight camera head: 105g

Accessories



300-001-602
High-End 1/3" MOS Camera Head
Full HD (included in the set)



300-001-615
Title generator keyboard, waterproof IP68
US-Layout, cable 1.8m, color: grey,
305x140x15mm, 600gr.



300-001-610
Footswitch for HD Cameras
cable length: 2m

TV COUPLER



4K
FULL HD

300-001-823 f=14mm
300-001-824 f=25mm
300-001-825 f=28mm
c-mount, soakable



4K
FULL HD

250-015-030
Zoom-TV-Adapter
f=15 - 30mm
c-mount, soakable



4K
FULL HD

250-015-036
Zoom-TV-Adapter
f=15 - 36mm
c-mount, soakable



FULL HD

300-001-813 f=14mm
300-001-814 f=25mm
300-001-815 f=28mm
c-mount, soakable



FULL HD

250-022-018 f=18mm
250-022-025 f=25mm
90°, c-mount, soakable

LED LIGHT SOURCE



Features:

- 6500K daylight color temperature
- 375 Watt equivalent output power
- > 20.000 hours lifetime
- Turret connector for the most common light cables (Olympus, Wolf, Gyrus ACMI, Storz)
- Manual or automatic light control
- Waterproof touch control panel
- Low heat conduction by improved LED technology
- Boost function: light increase of 20%

300-900-300 LUMEN LED 3 PRO Light Source, 110 - 230V



DELIVERY CONTAINS:

- Main Unit (300-900-300)
- BNC cable (optional)
- Mains cable
- Operating manual

Technical data

Line voltage	100-240 VAC
Video input	1x Video, BNC
Protection class	1 (acc. IEC 60601-1) 1 (acc. MDD 93/42/EEC)
Safety class	BF
Standard	IEC 60601-1-2
Dimensions	380 x 120 x 360mm (WxHxD)
Weight	approx. 11kg



A OLYMPUS

B RICHARD WOLF

C GYRUS ACMI

D STORZ
KARL STORZ ENDOSKOPE

Turret, multi connector


LED LIGHT SOURCE



Features:

- More than 20.000h lamp lifetime
- 180 Watt equivalent power output
- Manual or automatic light control
- Waterproof touch control panel
- Low heat conduction by improved LED technology
- Daylight color temperature 6000 – 6500° K for accurate color representation
- Boost function: light increase of 20%

300-900-180 LUMEN LED 1 Light Source




DELIVERY CONTAINS:

- Main Unit (300-900-180)
- BNC cable
- Mains cable
- Operating manual

Technical data

Line voltage	100-240 VAC
Power consumption	130VA
Video input	1x Video, 1x S-Video (Y/C) (automatic version only)
Protection class	1 (acc. IEC 60601-1) 1 (acc. MDD 93/42/EEC)
Safety class	BF
Dimensions	380 x 120 x 360mm (WxHxD)
Weight	approx. 11kg

HIGH TEMPERATURE LED/XENON LIGHT CABLES



Adaptor pieces to use with **endoscope**

XENON/LED Light Guide Cables with melt and bending protection

Adaptor pieces to use with **light source**

 250-042-000 RZ/Storz/Olympus	250-025-230 Ø 2.5 mm, 2300 mm 250-025-300 Ø 2.5 mm, 3000 mm 250-035-180 Ø 3.5 mm, 1800 mm 250-035-230 Ø 3.5 mm, 2300 mm 250-035-300 Ø 3.5 mm, 3000 mm 250-048-180 Ø 4.8 mm, 1800 mm 250-048-230 Ø 4.8 mm, 2300 mm 250-048-300 Ø 4.8 mm, 3000 mm	 250-039-000 RZ/Storz  250-038-000 ACMI/CUDA  250-048-000 Luxtec  250-051-000 Stryker  250-037-000 Wolf  250-044-000 Olympus  250-045-000 Olympus - Standard  250-052-000 Heine  250-059-000 Pentax
---	--	---



COMPATIBILITY NOTE:

A fiber optic cable with a large diameter should not be connected to a smaller diameter scope.
 For example: Ø 4.8mm fiber optic cables connected to a Ø 2.9mm may risk high temperature at the light entry of the endoscope and may risk patient safety.

ELECTROSURGICAL GENERATOR



Features:

- Monopolar/bipolar cutting and coagulation
- Automatic identification of connected instruments
- Powerstart System simplifies the beginning of cutting process
- Automatic power regulation
- 105 different memory locations in 7 groups
- Various output configuration of the sockets available
- 10" Touchscreen
- Tilttable display
- Footswitch (optional)
- Available with Argon and / or the RZeal
- Real-time output monitoring

DELIVERY CONTAINS:

- Main Unit
 - Mains Cable
- (Footswitch not included)

700-350-400 230V
700-350-410 110V
G400 ElectroSurgical Generator

Technical data

Line voltage	110-230 VAC
Power consumption	500W
Generator operation frequency	333kHz
Display	TFT 10"
Defibrillation impulse resistance	acc. to to EN 60601-1
NEM system	neutral electrode control optical indication 7 levels required resistance < 240 Ω
Gas type	pure argon 4.8 (99,998%) or higher
Inlet gas pressure	0.25-0.4 [MPa]
Monopolar coagulation / cutting	200W / 400W
Bipolar coagulation / cutting	400W / 400W

Gas outflow	0.1-10 l/min, regulated 0.1 l/min in whole range
Pressure range	reducer (0.4 MPa) on argon cylinder
Vessel sealing system	9 effects/300W
AutoStop	automatic stop of sealing process after achieving the optimal effect
Dimensions	495 x 415 x 225mm (LxWxH)
Weight	12.2kg

Accessories



700-350-405
2-Button-Footswitch, length: 5m
700-350-402
3-Button-Footswitch, length: 5m, quick access key



Disposable Neutral Electrode
hydrogel, split, for adults or children
10x5 pcs. per package
700-503-000 Adults
700-503-001 Children



Cable for Disposable Neutral Electrode
5m, RZ compatible only
700-503-003

ELECTROSURGICAL GENERATOR



Features:

- Monopolar/bipolar cutting and coagulation
- Automatic instrument detection system
- Intuitive Touchscreen
- Automatic power regulation
- Configurable working modes
- Real-time parameter control
- Special endoscopic modes
- Bipolar resection
- Easy software update

700-340-400 230V
700-340-410 110V
G44 Electrosurgical Generator

DELIVERY CONTAINS:

- Main Unit ■ Mains Cable
(Footswitch not included)

Technical data

Line voltage	110-240 VAC
Power consumption	1100VA
Generator operation frequency	333kHz
Display	TFT 7" WVGA
Spectrum result system	automatic power regulation
Monopolar coagulation	max. 200W (200Ω)
Monopolar cutting	max. 400W (100Ω)
Bipolar coagulation	max. 400W (25Ω)
Bipolar cutting	max. 330W (150Ω)

NEM system	neutral electrode control, optical indication, required resistance < 240 Ω
Dimensions	385 x 305 x 141mm (LxWxH)
Weight	7kg

Accessories



700-350-405
2-Button-Footswitch, length: 5m
700-350-402
3-Button-Footswitch, length: 5m, quick access key



Disposable Neutral Electrode
hydrogel, split, for adults or children
10x5 pcs. per package
700-503-000 Adults
700-503-001 Children



Cable for Disposable Neutral Electrode
5m, RZ compatible only
700-503-003

LITHOTRIPTER



Features:

- Pure mechanical energy transmission
- The smallest probe diameter of 0.6mm
- Auxiliary instruments can also be inserted
- Very simple handling and cleaning
- Compatible to other RZ and other instruments for the therapy of ureteric, renal and bladder stones
- Non-conductive hand piece to ensure safety for patients as well as users

DELIVERY CONTAINS:

- Main Unit (253-750-100)
- Operating manual ■ Mains cable

253-750-100 L100 Lithotripter



Technical data

Line voltage	115-230 VAC
Pneumatic power source	3.5 – 6.5 bar, compressed medical air
Pulse rate control (ppm)	50 – 600 ppm 10 ppm increments between 50 - 100 ppm 50 ppm increments between 100 - 600 ppm

Dimensions	254 x 113 x 305mm (WxHxD)
Weight	6.5kg

Accessories



Litho Probe, 605mm working length

253-850-006	0.6mm	253-850-016	1.6mm
253-850-008	0.8mm	253-850-020	2.0mm, 450mm length
253-850-010	1.0mm	253-850-025	2.5mm, 450mm length




253-750-150
Footswitch



253-750-200
Hand piece, anodised

L100 Lithotripter Modular System Chart




Handpiece + RZ Probes


Telescope + Sheaths

Auxiliary Instruments


URETEROSCOPY



253-750-200



253-850-006
/-008/-010/-012





253-905-540




253-954-005

CYSTOSCOPY




253-750-200



253-850-020/-025



253-900-101



253-900-111
253-900-116




253-900-210

PERCUTANEOUS NEPHROSCOPY




253-750-200



253-850-020/-025



253-900-100



253-900-110
253-900-115



253-900-205



300-001-168
32" HD Monitor

Technical data

Screen size	32" diagonal TFT LCD
Resolution	1920 x 1080
Aspect ratio	16:9
Viewing angle	178°
Numbers of colors	16.7 Million (8-bit)
Video Inputs	2x HDMI, 1x DP, 2x SDI
Video Outputs	1x HDMI, 2x SDI
Weight	11.5kg



300-001-232
Radiance 32" HD Monitor

300-010-226
Radiance 27" HD Monitor

Technical data


Screen size	32" / 27" diagonal
Resolution	1920 x 1080
Aspect ratio	16:9
Viewing angle	178°
Numbers of colors	1.07 Billion (10-bit)
Video Inputs	DVI-I and 3G-SDI Input/Output
Video Outputs	DVI-D
Weight	11.0kg (300-001-232) 8.9kg (300-010-226)



300-001-031
UHD 4K Monitor, 31"

Technical data

Screen size	31" diagonal
Resolution	4096 x 2160
Aspect ratio	17:9
Viewing angle	178°
Number of colors	1.08 Billion (10-bit)
Video Inputs	1 x HDMI 2.0 1 x DP 1.2 (SST) & 1 x DP 1.2 (MST) 1 x DVI (single link), 4 x SDI (3G)
Video Outputs	1 x DP 1.2 (SST) 1 x DVI (single link), 4 x SDI (3G)
Weight	10.5kg



300-001-055
UHD 4K Monitor, 55"

Technical data

Screen size	55" diagonal
Resolution	3840 x 2160
Aspect ratio	16:9
Viewing angle	178°
Numbers of colors	1.08 Billion (10-bit)
Video Inputs	1 x HDMI 2.0, 2 x DP 1.2 (55T) 1 x DVI (single link), 2 x DVI (dual link) 4 x SOi (3G), 2 x SOi (12G)
Video Outputs	1 x DVI (single link), 1 x DP 1.2 (SST) 4 x SDI (3G), 2x SDI (12G)
Weight	37.35kg



300-001-028
Monitor Table Stand

Technical data




Column height	305mm, 12"
Oriantations	Landscape or page mode
Adjustments	+/- 20°
Max. load	31" / 13.6kg / 30lbs

All monitors will be delivered without table stand - need to be ordered separately.


RESECTION SIMULATOR / TRAINER

Features:

- Essential hand movements can be learned
- Anatomical substrate can be detected and removed by using the resectoscope
- Usage in workshops, exhibitions or at the customer's location



900-253-000 GIUSEPPE Urological Simulator / Trainer, prostate and bladder, incl. pump and accessories, without substrates



Urological Substrate with and without anatomical structures

- Realistic tissue
- Myomas and Polyps structures
- Close to reality environment

	Quantity	Storable
900-253-001 Anatomical structure of the prostate	12	6 months
900-253-002 Anatomical structure of the bladder	6	3 months
900-253-003 Without anatomical structure	12	12 months

250-015-030	136	253-000-219	25	253-000-545	8, 59
250-015-036	136	253-000-220	25	253-000-547	59
250-020-025	136	253-000-222	25	253-000-555	58
250-020-030	136	253-000-225	25	253-000-666	54, 71
250-025-230	139	253-000-250	54	253-000-676	54
250-025-300	139	253-000-251	54, 71	253-000-677	54
250-035-180	139	253-000-256	5, 7, 12, 47, 68,	253-000-678	54
250-035-230	139		69, 70, 74, 77, 86, 87	253-000-724	47
250-035-300	139	253-000-258	86, 87	253-000-850	117
250-037-000	139	253-000-266	47, 69	253-000-900	117
250-038-000	139	253-000-300	68, 69	253-000-901	117
250-039-000	139	253-000-301	10, 11, 69, 70,	253-000-950	2, 5, 13, 117
250-040-000	139		72, 73	253-000-951	117
250-041-000	139	253-000-302	12, 69, 70, 72, 73	253-000-975	117
250-042-000	139	253-000-305	68, 70	253-000-976	117
250-044-000	139	253-000-312	12, 68, 69	253-200-007	104
250-045-000	139	253-000-316	68, 72	253-200-009	104
250-048-000	139	253-000-318	11, 68, 73	253-200-011	105
250-048-180	139	253-000-319	10, 68, 72	253-200-075	3, 105
250-048-230	139	253-000-320	70	253-200-095	3, 105
250-048-300	139	253-000-321	70	253-200-100	110
250-051-000	139	253-000-322	70	253-200-128	111
250-052-000	139	253-000-323	70	253-200-200	4, 111
250-059-000	139	253-000-330	70	253-200-210	4, 111
253-000-044	48	253-000-331	70	253-200-228	111
253-000-054	48	253-000-332	70	253-200-328	111
253-000-082	72	253-000-333	70	253-200-428	111
253-000-084	5, 48	253-000-350	71	253-200-500	110
253-000-094	7, 48	253-000-352	8, 10, 58, 68, 72	253-200-600	4, 110
253-000-100	26	253-000-356	47, 74, 86, 87	253-202-003	3, 104, 106
253-000-101	26	253-000-358	86, 87	253-203-128	4, 111
253-000-102	2, 26	253-000-362	72	253-203-228	111
253-000-105	2, 26	253-000-366	47	253-203-328	4, 111
253-000-106	26	253-000-400	7, 46	253-203-428	4, 111
253-000-107	26	253-000-401	46	253-204-003	3, 104, 106
253-000-108	26	253-000-410	8, 58	253-204-128	111
253-000-117	25	253-000-411	58	253-204-228	111
253-000-119	25	253-000-501	59	253-204-328	111
253-000-120	25	253-000-502	59	253-204-428	111
253-000-122	25	253-000-503	59	253-219-100	6, 14, 105, 108,
253-000-125	25	253-000-505	59		109
253-000-200	5, 13, 45, 79	253-000-507	59	253-219-130	3, 105, 108, 109
253-000-201	45	253-000-524	47	253-220-011	106
253-000-210	57	253-000-541	8, 59	253-220-014	106
253-000-211	57	253-000-542	59	253-220-025	106
253-000-217	25	253-000-543	8, 59	253-220-030	106

253-221-000	104	253-300-222	49	253-334-335	120
253-221-030	104	253-300-251	55	253-335-225	121
253-227-000	106	253-300-252	55	253-335-335	120
253-227-030	106	253-300-253	55	253-336-225	121
253-227-070	106	253-300-255	55	253-336-335	120
253-250-200	6, 14, 108, 109	253-300-256	55	253-337-225	121
253-250-211	6, 108	253-300-258	55	253-337-335	120
253-250-213	108	253-300-259	55	253-340-335	120
253-251-111	6, 108	253-300-999	5, 7, 8, 9, 13, 80	253-341-335	120
253-251-112	6, 108	253-310-014	107	253-342-335	120
253-251-113	6, 108	253-310-015	107	253-343-335	120
253-251-114	108	253-310-100	107	253-344-335	120
253-251-115	6, 108	253-310-101	107	253-345-335	120
253-252-111	108	253-310-107	107	253-346-335	120
253-252-112	108	253-310-390	119	253-347-335	120
253-252-113	108	253-310-392	119	253-348-335	120
253-252-114	108	253-315-000	25	253-349-335	120
253-252-115	108	253-320-225	121	253-350-335	120
253-255-001	14, 109	253-320-335	120	253-351-335	120
253-265-010	14, 109	253-321-225	121	253-352-335	120
253-265-101	14, 109	253-321-335	120	253-353-335	120
253-265-102	14, 109	253-322-225	121	253-354-335	120
253-265-103	109	253-322-335	120	253-355-335	120
253-265-104	14, 109	253-323-225	121	253-356-335	120
253-300-017	2, 24	253-323-335	120	253-357-335	120
253-300-019	24	253-324-225	121	253-370-003	37, 111
253-300-020	24	253-324-335	120	253-370-004	111
253-300-022	2, 24	253-325-225	121	253-370-005	2, 37, 111
253-300-025	19, 24, 86	253-325-335	120	253-370-007	2, 37
253-300-201	5, 49, 51	253-326-225	121	253-370-009	121
253-300-202	49	253-326-335	120	253-370-010	37
253-300-203	5, 49, 51	253-327-225	121	253-370-021	121
253-300-204	49	253-327-335	120	253-370-103	37
253-300-205	5, 49, 51	253-328-225	121	253-370-253	3, 4, 104, 106
253-300-206	5, 49, 51	253-328-335	120	253-371-003	37, 111
253-300-207	49	253-329-225	121	253-371-005	37, 111
253-300-208	49	253-329-335	120	253-371-007	37
253-300-209	50	253-330-225	121	253-371-253	3, 4, 104, 106
253-300-210	50	253-330-335	120	253-372-005	37
253-300-211	50	253-331-225	121	253-372-007	37
253-300-212	50	253-331-335	120	253-375-015	121
253-300-213	50	253-332-225	121	253-400-200	52
253-300-214	50	253-332-335	120	253-400-201	7, 52
253-300-215	50	253-333-225	121	253-400-202	52
253-300-218	50	253-333-335	120	253-400-203	7, 52
253-300-219	50	253-334-225	121	253-400-204	52

253-400-205	7, 52	253-500-315	81	253-755-500	85
253-400-206	52	253-500-321	13, 79	253-760-000	19, 86
253-400-211	52	253-500-325	13, 79	253-780-000	87
253-400-222	52	253-500-421	79	253-800-000	114
253-400-501	64	253-500-500	81	253-800-900	115
253-400-502	64	253-500-561	13, 80	253-804-000	13, 19, 23, 44,
253-400-503	64	253-500-562	80		68, 79, 84
253-400-504	64	253-500-563	13, 80	253-804-012	23, 44, 68, 69,
253-400-505	64	253-500-564	13, 80		70, 76, 77
253-400-506	64	253-500-581	81	253-804-030	2, 5, 7, 12, 19,
253-400-508	64	253-500-582	81		23, 44, 68, 69, 70, 76, 77, 85, 86, 87
253-400-601	63, 112	253-580-007	36	253-804-070	23, 85
253-400-602	63, 112	253-581-007	36	253-804-400	22
253-400-606	63, 112	253-581-045	35	253-804-412	22
253-400-701	53	253-582-045	35	253-804-430	22
253-400-702	53	253-582-055	35	253-804-470	22
253-400-703	53	253-582-065	35	253-804-530	54, 71
253-400-704	53	253-584-005	36	253-829-000	81, 107
253-420-100	2, 34	253-585-005	2, 15, 35, 102	253-829-012	72, 73
253-421-100	34	253-610-001	37, 63, 112	253-829-030	9, 72, 73, 107
253-422-100	34	253-610-002	37, 63, 112	253-850-006	20, 142, 143
253-423-100	34	253-610-003	2, 3, 5, 6, 37, 63,	253-850-008	20, 142
253-425-100	34		112	253-850-010	20, 142
253-426-100	34	253-610-004	37, 63, 112	253-850-012	20
253-427-100	34	253-650-010	120	253-850-016	20, 142
253-428-100	12, 74	253-650-012	120	253-850-020	20, 142, 143
253-428-500	47	253-650-014	120	253-850-025	142
253-429-100	34	253-650-016	120	253-870-023	76
253-430-100	86	253-650-018	120	253-872-000	76
253-431-100	34	253-650-020	120	253-873-000	76
253-500-000	2, 5, 19, 117	253-650-022	120	253-874-000	69, 77
253-500-001	117	253-650-024	120	253-900-000	92
253-500-003	117	253-650-026	120	253-900-010	89, 102
253-500-004	117	253-650-028	120	253-900-100	93, 94, 143
253-500-010	117	253-650-030	120	253-900-101	93, 94, 143
253-500-020	117	253-650-032	120	253-900-105	93, 94
253-500-100	13, 82	253-650-034	120	253-900-110	93, 143
253-500-101	82	253-650-036	120	253-900-111	93, 143
253-500-102	82	253-750-010	142, 143	253-900-115	93, 143
253-500-103	82	253-750-100	20, 142, 143	253-900-116	93, 143
253-500-108	82	253-750-150	20, 142, 143	253-900-120	94
253-500-115	81	253-750-200	20, 142, 143	253-900-121	12, 18, 68, 74, 94
253-500-125	81	253-755-010	19, 84	253-900-130	18, 94
253-500-221	13, 79	253-755-020	19, 84	253-900-131	94
253-500-229	79	253-755-030	19, 84	253-900-135	18, 94
253-500-230	13, 79	253-755-040	84	253-900-136	94

253-900-140	16, 89	253-900-921	9	253-952-505	31, 100
253-900-141	16, 90	253-900-922	14	253-953-003	31, 100
253-900-142	16, 90	253-900-923	17	253-953-004	31, 100
253-900-143	16, 90	253-900-924	18	253-953-005	31, 100
253-900-145	16, 90	253-900-925	19	253-953-503	31, 100
253-900-146	16, 90	253-900-926	20	253-953-504	31, 100
253-900-147	16, 90	253-901-220	11, 68, 75	253-953-505	31, 100
253-900-148	16, 90	253-904-220	10, 11, 75	253-954-003	31, 101
253-900-150	93	253-905-000	102, 125	253-954-004	31, 101
253-900-151	93	253-905-220	10, 11, 68, 75	253-954-005	15, 31, 101, 143
253-900-155	18, 95	253-905-400	89, 102	253-954-455	91
253-900-156	18, 95	253-905-530	99	253-954-503	31, 101
253-900-157	18, 95	253-905-540	15, 99, 143	253-954-504	31, 101
253-900-160	94	253-910-000	132	253-954-505	31, 101
253-900-161	94	253-910-001	132	253-955-003	31, 101
253-900-171	90	253-910-002	132	253-955-004	31, 101
253-900-173	90	253-910-003	132	253-955-005	31, 101
253-900-175	90	253-910-530	99	253-955-503	31, 101
253-900-200	91, 95	253-910-540	99	253-955-504	31, 101
253-900-201	18, 95	253-915-530	99	253-955-505	31, 101
253-900-202	18, 95	253-915-540	99	253-956-003	31, 101
253-900-205	96, 143	253-920-090	121	253-956-004	31, 101
253-900-206	96	253-920-105	121	253-956-005	15, 31, 101
253-900-210	96, 143	253-920-120	121	253-956-455	16, 17, 91
253-900-211	96	253-920-135	121	253-956-503	31, 101
253-900-212	96	253-920-150	121	253-956-504	31, 101
253-900-213	96	253-940-005	35, 102	253-956-505	31, 101
253-900-245	96	253-950-003	30, 100	253-957-003	31, 101
253-900-246	96	253-950-004	30, 100	253-957-004	31, 101
253-900-250	18, 96	253-950-005	15, 30, 100	253-957-005	31, 101
253-900-251	96	253-950-455	16, 17, 91	253-957-105	27
253-900-252	18, 96	253-950-503	30, 100	253-957-503	31, 101
253-900-253	96	253-950-504	30, 100	253-957-504	31, 101
253-900-815	92	253-950-505	30, 100	253-957-505	31, 101
253-900-830	17, 92	253-951-003	30, 100	253-958-003	32, 101
253-900-888	17, 92	253-951-004	30, 100	253-958-004	32, 101
253-900-900	2	253-951-005	30, 100	253-958-005	32, 101
253-900-901	5	253-951-503	30, 100	253-958-105	27
253-900-902	13	253-951-504	30, 100	253-958-455	91
253-900-903	15	253-951-505	30, 100	253-958-503	32, 101
253-900-905	6	253-952-003	31, 100	253-958-504	32, 101
253-900-909	16	253-952-004	31, 100	253-958-505	32, 101
253-900-911	7	253-952-005	15, 31, 100	253-959-003	32, 101
253-900-918	3	253-952-455	16, 17, 91	253-959-004	32, 101
253-900-919	4	253-952-503	31, 100	253-959-005	32, 101
253-900-920	12	253-952-504	31, 100	253-959-503	32, 101

253-959-504	32, 101	301-701-030	98, 118	351-100-500	30, 100
253-959-505	32, 101	301-701-060	98, 118	351-100-550	32, 100
253-960-025	102	301-701-100	38, 118	351-107-107	27
253-960-360	91	320-603-455	119	351-107-109	27
253-961-102	38	320-720-080	119	351-114-105	27
253-961-103	38	320-720-100	119	351-114-107	2, 27
253-962-003	102	320-720-120	119	351-114-109	27
253-990-800	38	320-720-140	119	351-120-105	27, 80
300-001-028	145	320-720-160	119	351-120-107	2, 27
300-001-031	145	320-720-180	119	351-120-109	27
300-001-055	145	320-720-200	119	351-124-107	27
300-001-158	144	320-720-220	119	351-400-200	9, 61
300-001-168	144	320-720-240	119	351-400-201	61
300-001-232	144	320-721-080	119	351-829-000	56, 60, 68
300-001-600	135	320-721-100	119	351-829-012	56, 60, 68
300-001-601	135	320-721-120	119	351-829-030	8, 10, 11, 56, 60,
300-001-602	135	320-721-140	119		68
300-001-610	133, 134, 135	320-721-160	119	351-905-105	118
300-001-615	133, 134, 135	320-721-180	119	600-050-000	130
300-001-800	134	320-721-200	119	600-100-000	130
300-001-802	134	320-721-220	119	600-255-500	129
300-001-804	134	320-721-240	119	600-255-530	129
300-001-813	136	320-723-080	119	600-255-535	129
300-001-814	136	320-723-100	119	600-255-550	129
300-001-815	136	320-723-120	119	600-255-555	129
300-001-823	136	320-723-140	119	600-255-565	129
300-001-824	136	320-723-160	119	600-255-575	129
300-001-825	136	320-723-180	119	600-255-600	129
300-001-840	133	320-723-200	119	600-255-605	129
300-001-841	133	320-723-220	119	600-281-001	129
300-010-226	144	320-723-240	119	600-281-002	128
300-011-178	118	320-998-000	121	600-281-003	127
300-011-179	118	320-999-000	121	600-281-004	127
300-011-181	118	351-000-155	61	600-292-052	125
300-011-184	10, 11, 12, 69,	351-000-180	9, 61	600-405-500	128
	70, 72, 73, 118	351-000-185	9, 11, 61, 68, 73	600-405-530	128
300-011-506	118	351-000-186	61	600-405-535	128
300-011-515	118	351-000-195	61, 73	600-405-550	128
300-135-905	130	351-000-201	9, 62	600-405-555	128
300-135-910	130	351-000-202	62	600-405-570	128
300-215-000	38, 130	351-000-203	9, 62	600-405-575	128
300-400-900	118	351-000-205	9, 62	600-405-600	128
300-400-901	118	351-000-206	62	600-405-605	128
300-420-065	130	351-100-105	2, 27	600-450-070	125
300-900-180	138	351-100-107	27	600-460-055	125
300-900-300	137	351-100-109	27	600-460-056	125

600-462-052	125	700-900-220	130
600-464-229	124	700-900-250	130
600-485-500	127	900-253-000	146
600-485-530	127	900-253-001	146
600-485-535	127	900-253-002	146
600-485-550	127	900-253-003	146
600-485-555	127		
600-485-570	127		
600-485-575	127		
600-485-600	127		
600-485-605	127		
600-501-000	123		
600-501-002	123		
600-505-000	123		
600-505-002	123		
600-506-002	123		
600-508-000	124		
600-511-000	124		
600-540-500	127		
600-540-530	127		
600-540-535	127		
600-540-550	127		
600-540-555	127		
600-540-570	127		
600-540-575	127		
600-540-600	127		
600-540-605	127		
600-600-191	102		
600-605-220	126		
600-605-411	126		
600-662-077	102, 125		
600-662-540	102		
600-950-100	130		
660-166-100	130		
700-340-400	141		
700-340-410	141		
700-350-400	140		
700-350-402	140, 141		
700-350-405	140, 141		
700-350-410	140		
700-500-453	7, 8, 9, 63, 112		
700-501-453	37, 63, 112		
700-503-000	140, 141		
700-503-001	140, 141		
700-503-003	140, 141		
700-900-160	130		

