


SURFACE OPTIONS



RZ
Medizintechnik
GmbH


Progress in Surgery

Matt finish



Brushed, slightly glossy surface.
Standard finish.

Sand-blasted



Matt, smooth surface.
Optional finish.

Mirror finish



Highly glossy and smooth surface.
Optional finish.

SCISSORS
IDENTIFICATION



CC
CERACUT®
Overall Ceramic coating



TC
Tungsten Carbide Inserts
Golden ring, matt finish



SC
Super Cut
Black coated ring,
Mirror finish



TC CC
Tungsten Carbide Inserts
+ Ceramic coated surface
Golden ring



TC SC
Tungsten Carbide Inserts
Super Cut
Golden ring,
Sand blasted surface