



Laryngo-Pharyngoscope Retractor



Designed for precision and control.

OSTRA Laryngo-Pharyngoscope Retractor



Proven concept – now with RZ Quality

The new OSTRA retractor system from RZ Medizintechnik brings a trusted and established design back into the operating room. With precision engineering and consistent quality, OSTRA ensures that surgeons can rely on a complete and dependable solution for both robotic and conventional transoral surgery.

OSTRA integrates into a comprehensive ENT portfolio, offering not only a single product, but an entire system of coordinated instruments for safe and efficient airway and laryngo-pharyngeal procedures.

Key features and advantages

- **Wide range of blades**
Tongue, vallecular, and laryngeal
- **Adaptable design**
Accommodates different anatomical conditions and surgical techniques
- **Open lateral frame**
Enables free movement for robotic instruments and camera, ensuring an unobstructed view
- **Independent blade adjustment**
Precise control of angle and depth for customized exposure.
- **Secure fixation**
Fully compatible with the RZ Retractor Holder Set



RZ Medizintechnik GmbH

Unter Hasslen 20/22
78532 Tuttlingen, Germany
Phone: +49 7462 9470-0
sales@rz-medizintechnik.com
www.rz-medizintechnik.com

430-326-000

OSTRA Laryngo-Pharyngoscope Retractor
Complete Set

