




*progress in surgery*

# ArthVANCE<sup>®</sup>

## Suture Passer





# ARTHVANCE® SUTURE PASSER


## Rotator Cuff Repair

Simultaneously one-handed passing and retrieving of suture

- Efficient tissue bite (wide jaw opening and deep bite)
- Minimal perforation of tissue (small hole)
- Reusable instrument and reliable construction




Suture Passer with trap door frontal opening **without ratchet**  
170mm working length 225-820-217




Suture Passer with trap door frontal opening **with ratchet**  
170mm working length 225-820-213



Needle with lateral opening  
0.3 x 240mm 225-820-219



Automatic trap loop hold




Slim design



Ratchet mechanism (225-820-213)



Strong grip 4-finger trigger



Luer lock flush port (225-820-213)



*progress in surgery*

## RZ Medizintechnik GmbH

Unter Hasslen 20  
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50

[sales@rz-medizintechnik.com](mailto:sales@rz-medizintechnik.com)  
[www.rz-medizintechnik.com](http://www.rz-medizintechnik.com)

