

# Laparoscopy Hand Instruments





# LAPAROSCOPY HAND INSTRUMENTS, 360° ROTATABLE



## THE STANDARD

The modular system for whatever you need

- Long-life and reliable materials Large variety of handle designs Ratchet mechanism / deactivation

More than 400 different laparoscopic inserts (Ø 3.0 / 3.5 / 5.0 / 10mm)



The RZ spherical head locking system for fast assembling / disassembling



Modular 3-piece detachable instrument system consists of:



#### Please find more information in our Laparoscopy catalogue





Outer Sheath

Insert

# LAPAROSCOPY HAND INSTRUMENTS, 360° ROTATABLE



# THE ENDO TIP SYSTEM

The environmental solution for the always sharp scissors

- Reusable handle and shaft, disposable tip Always having a sharp tip
- Ecological and economical reduces waste and saves budget



Complete Instrument consists of:

**D** PEEK Polymer Handle (300-450-300)

**2** Actuator rod (300-450-006)

3 Outer tube including screw cap blue, Ø 5mm (300-450-005)

Insert

#### PEEK POLYMER HANDLES





The EndoTip system allows to reuse the handle and the sheath of the Instrument, but gives the opportunity to easily change the tip. That provides an always sharp scissors without creating the waste of a complete single use instrument.

#### 4



**METZENBAUM** non-sterile, limited reusable 300-450-177 medium, 15mm curved 300-450-171 short, 10mm, curved 300-450-173 long, 18 mm, curved 300-450-174 long, 18 mm, straight 300-450-176 medium, 15mm, straight



**METZENBAUM** sterile, single use 300-450-177S medium, 15mm curved



# LAPAROSCOPY HAND INSTRUMENTS, 360° ROTATABLE



# SINGLE USE SHEATH SYSTEM

Convenient and sharp

- Combine your preferred handle with a sharp singe use scissors No need for resharpening of your most preferred Metzenbaum High quality handle combined with single use scissors •



The fully disposable sheath system gives you the opportunity to combine your prefered MIC instrument handle with a disposable sheath. That safes space on your tray and reduces time in reprocessing of the equipment.

#### METZENBAUM

Scissors Inserts, sterile, box with 10 pcs., Ø 5mm, 330mm working length medium, 15mm curved

300-600-177S

METZENBAUM

non-sterile, limited reusable, 1 pc.	
straight, short, 10mm	300-600-17
curved, short, 10mm	300-600-17
curved, 18 mm, long	300-600-17
straight, 18 mm, long	300-600-17
straight, 15mm, medium	300-600-17
curved, 15mm, medium	300-600-17

PEEK Polymer Handle, with HF connection 300-750-101 without ratchet





### RZ Medizintechnik GmbH

Unter Hasslen 20/22 78532 Tuttlingen, Germany

Phone: +49 7462 9470-0

sales@rz-medizintechnik.com www.rz-medizintechnik.com

