




## Laryngo-Pharyngoscope Retractor



Designed for precision and control.

# OSTRA Laryngo-Pharyngoscope Retractor



## Proven concept – now with RZ Quality

The new OSTRa retractor system from RZ Medizintechnik brings a trusted and established design back into the operating room. With precision engineering and consistent quality, OSTRa ensures that surgeons can rely on a complete and dependable solution for both robotic and conventional transoral surgery.

OSTRa integrates into a comprehensive ENT portfolio, offering not only a single product, but an entire system of coordinated instruments for safe and efficient airway and laryngo-pharyngeal procedures.

## Key features and advantages

- **Wide range of blades**  
Tongue, vallecular, and laryngeal
- **Adaptable design**  
Accommodates different anatomical conditions and surgical techniques
- **Open lateral frame**  
Enables free movement for robotic instruments and camera, ensuring an unobstructed view
- **Independent blade adjustment**  
Precise control of angle and depth for customized exposure.
- **Secure fixation**  
Fully compatible with the RZ Retractor Holder Set



### RZ Medizintechnik GmbH

Unter Hasslen 20/22  
78532 Tuttlingen, Germany  
Phone: +49 7462 9470-0  
sales@rz-medizintechnik.com  
www.rz-medizintechnik.com

430-326-000

OSTRa Laryngo-Pharyngoscope Retractor  
Complete Set

